



## EM - TECHCOLOR PIPETTES

Pipettes are used for accurate measurements and transfer of liquids. Volumetric pipettes will only deliver specific quantities. Graduated pipettes allow take - up of different liquid volumes and delivery of the same or different quantities. Graduated and volumetric pipettes are of EM - TECHCOLOR brand. They are calibrated to deliver ("Ex") at the reference temperature of 20°C. In the case of calibrated pipettes and Class 'AS' pipettes, the tolerances on volume content corresponds to the accuracy limits of **ISO** and **DIN** standards. The accuracy limits for EM - TECHCOLOR Class 'B' type are about one

and half times that of Class 'AS'. This means that they are better than required by **ISO** and **DIN**. By specifying a Class 'AS' within the frame work or their 15th modifications order, the Federal German weights and measures authorities have recognised that an overwhelming portion of volumetric measurements, particularly in clinical laboratories is carried out on water or on dilute aqueous solutions. Thus, apparatus with the same accuracy limits, but with run-out time considerably shorter than hitherto demanded, are now admitted for calibration.

### Tolerances for Graduated measuring Pipettes calibrated to deliver

Capacity ml	Accuracy Limits Class 'AS' suitable for official calibration		Accuracy Limits Class 'B'	
	ISO 835 DIN 12697 ± ml	ISO 835 DIN 12695 ± ml	EM-TECHCOLOR ± ml	
0.1 <sup>1</sup>	-	-	0.0015	
0.2 <sup>1</sup>	-	-	0.0030	
1	0.007	0.01	0.0100	
2	0.010	0.02	0.0150	
5	0.030	0.05	0.0450	
10	0.050	0.10	0.0750	
25	0.100	0.20	0.1500	

1 Non - DIN, Non - ISO size, Graduated Pipettes 0.1 and 0.2 ml are calibrated to contain ("In")

### Tolerances for Volumetric Pipettes calibrated to deliver

Capacity ml	Accuracy Limits Class 'AS' suitable for official calibration		Accuracy Limits Class 'B'	
	DIN 12691 ± ml	ISO 648 DIN 12690 ± ml	EM-TECHCOLOR ± ml	
5	0.015	0.03	0.0230	
10	0.020	0.04	0.0300	
20	0.030	0.06	0.0450	
25	0.030	0.06	0.0450	
50	0.050	0.10	0.0750	


**Pipettes, Graduated, Clear Glass, Short - Line, Graduated to Tip,  
Amber Stain Graduation, (Serological type), Class 'B', RIVIERA Brand**

Cat	Capacity	Graduation	Tolerance	Total	Qty. Per
No.	ml	Division ml	± ml	Length ±5mm	Pack Pcs
72600 14	0.10	0.01	0.01	320	10
72600 21	0.20	0.02	0.01	320	10
72600 29	0.50	0.05	0.01	320	10
72600 33	1.00	0.01	0.01	350	10
72600 41	2.00	0.02	0.02	350	10
72600 53	5.00	0.05	0.05	350	10
72600 60	10.0	0.10	0.10	350	10
72600 70	25.0	0.20	0.20	470	10



# VOLUMETRIC GLASSWARE FROM HIRSCHMANN LABORGERATE

GRADUATED PIPETTE, CLASS 'B'  
EM-TECHCOLOR




In any laboratory volume measurement is a continuous process. The products that aid such volume measurements, are Burettes, Graduated Pipettes, Volumetric Pipettes, Volumetric Flasks, and Measuring Cylinders. These aids are available in both class A and Class B and the selection of the class depends on the application and required accuracy.

All the volumetric products are either to contain or to deliver the liquids. To contain is signified by the mark shown on the ware "in" and to deliver is signified by the mark on the ware "Ex". Products such as Measuring Cylinders and Volumetric Flasks are meant for containing whereas products such as Pipettes and Burettes are meant for delivering. The Capacity of the "in" products will always be equivalent to the capacity printed on the ware. However the "Ex" products will show the capacity printed as that which is equivalent to the volume delivered. The residue adhering to the wetted glass surface is already taken into account while calibrating.

The division of volumetric glassware into classes A, AS, and B indicates the suitability for certification, degrees of accuracy and DELIVERY and WAITING TIMES. DELIVERY TIME is the time required for the liquid meniscus to pass from the calibration mark to standstill at the tip or at the second calibration mark. WAITING TIME starts when the liquid meniscus in the


VOLUMETRIC PIPETTE, CLASS 'AS'  
EM-TECHCOLOR




tip or at the second calibration mark has come to a stop. During the waiting time, residual liquid flows downward from the glass wall. This causes a new rise of the meniscus. After the waiting time, the pipette tip is wiped on the wall of the vessel. A small amount of liquid remains in the narrow section of the tip. But as explained previously this has been taken into account during calibration. Products of Class A and Class AS are Conformity Certified wares. Those, which are supposed to contain ("in"), will have class A and those, which are supposed to deliver ("Ex"), will have Class AS. Class AS has all the features of Class A and in addition its delivery is swift (S stands for swift delivery). Class AS will have a delivery time of about 15 seconds and a waiting period of another 15 seconds. The tolerances for accuracy for Class B remain within twice the limits defined by DIN and ISO for Class A/AS.

## The meniscus - How does it arise, how is it read ?


If the liquid molecules are attracted more strongly by the molecules of the glass wall (adhesion) than by their counterparts (cohesion), a concave meniscus is formed. The edge of the liquid level is drawn slightly upwards. If the diameter of a pipette is narrow enough (e. g. in a capillary), the adhesion force is sufficient to pull up not just the outer edge, but the entire liquid level. This is the whole "secret" of the capillary effect. If however, the cohesion force of a liquid is greater than the adhesion force of the glass wall, a convex meniscus is formed (e. g. Mercury)



The Schellbach stripe is a slender blue stripe in the middle of a white milkglass stripe printed at the back of the volumetric glassware. The light refraction on the surface of the liquid causes a narrowing (constriction) of the blue stripe between the top and bottom meniscus. According to DIN 12700, read-off is at the point of contact of the two tips



In the case of a convex meniscus, e. g. mercury, the German Industrial Standard DIN/ISO 384 and DIN 12700 requires the volume to be read at the highest point of the liquid's surface. The highest point of the meniscus must cover the lower edge of the graduation mark.



In the case of a concave meniscus, e. g. water, the German Industrial Standard DIN/ISO 384 and DIN 12700 requires the volume to be read at the lowest point of the liquid's surface. The lowest point of the meniscus must cover the upper edge of the graduation mark.

## Recognition of DAkkS under Indian Regulation of NABL


Hirschmann Laborgerate GmbH & Co. KG has been given Accreditation by DAkkS (Deutsche Akkreditierungsstelle GmbH) as competent to carry out calibrations. This accreditation implies that Hirschmann has the right to issue internationally recognized certificates for volumetric measuring instruments. Therefore, the tested quality of each product is not only documented but the certification represents the rigorous testing procedures as well as the excellence of the quality management systems. The DAkkS calibration certificate includes all measured values and measuring uncertainty of the tested device based on international norms. DAkkS is signatory to the Multilateral Agreements of EA, ILAC and IAF for Mutual Recognition.

DAkkS and NABL are connected and mutually recognized through the ILAC Organization. ILAC (International Laboratory Accreditation Cooperation) is the international organization for accreditation bodies operating in accordance with ISO/IEC 17011 and involved in the accreditation of Conformity Assessment Bodies including Calibration Laboratories (using ISO/IEC 17025), Testing Laboratories (using ISO/IEC 17025), Medical Testing Laboratories (using ISO 15189) and Inspection Bodies (using ISO/IEC 17020).

Accreditation allows you to make an informed decision when selecting an organization to carry out your testing, calibration, or inspection activities, as it demonstrates competence, impartiality and capability of that organization. As a result, accreditation helps to underpin the credibility, safety and performance of goods and services.

Accreditation bodies around the world, which have been evaluated by peers as competent, have signed an arrangement that enhances the acceptance of products and services across national borders. The purpose of this arrangement, the ILAC Mutual Recognition

Arrangement (MRA), (often referred to as ILAC Arrangement) is to create an international framework to support international trade through the removal of technical barriers.



In many economies there is an accreditation body recognized by government and industry (in India NABL) to carry out the assessment and verification against international standards of testing, calibration, inspection and certification activities in both the private and public sectors. ILAC is an organization that counts as its members laboratory accreditation bodies, regional organizations and stakeholders representing over 110 economies. The ILAC MRA allows making use of a global network of testing and calibration laboratories and inspection bodies that have been accredited to provide accurate and reliable results. The ILAC MRA acts as an internationally recognized 'stamp of approval' to demonstrate compliance against agreed standards and requirements. Accreditation and the ILAC MRA help regulators meet their own legislated responsibilities by providing a globally recognized system to accept accredited test and inspection reports.

### Signatories to the ILAC Mutual Recognition Arrangement

No.	Accreditation Body	Economy	Scope	Original Signing Date
29	Comite Francais d'Accreditation (COFRAC)	France	38 Testing ISO/IEC 17025 & ISO 15189 Calibration ISO/IEC 17025 Inspection ISO/IEC 17020	02 Nov 2000 02 Nov 2000 24 Oct 2012
30	<sup>(6)</sup> Deutsche Akkreditierungsstelle GmbH (DakKS)	Germany	39 Testing ISO/IEC 17025 & ISO 15189 Calibration ISO/IEC 17025 Inspection ISO/IEC 17020	02 Nov 2000 02 Nov 2000 24 Oct 2012
31	<sup>(7)</sup> Hellenic Accreditation System S.A. (ESYD)	Greece	40 Testing ISO/IEC 17025 & ISO 15189 Calibration ISO/IEC 17025	22 May 2004 22 May 2004
34	National Accreditation Board for Testing and Calibration Laboratories (NABL)	India	43 Testing ISO/IEC 17025 & ISO 15189 Calibration ISO/IEC 17025	02 Nov 2000 02 Nov 2000
35	National Accreditation Board for Certification Bodies (NABCB)	India	43 Inspection ISO/IEC 17020	16 Sept 2013
36	National Accreditation Body of Indonesia (KAN)	Indonesia	44 Testing ISO/IEC 17025 Testing ISO 15189 Calibration ISO/IEC 17025	20 June 2001 14 Mar 2013 30 Dec 2003



**Pipettes, Graduated, Clear Glass, Short - Line, Graduated to Tip, Amber Stain Graduation, (Serological type), Class 'B', EM-TECHCOLOR Brand** DIN EN ISO 835

Cat No.	Capacity ml	Graduation Division ml	Tolerance $\pm$ ml	Total Length $\pm$ 5mm	Qty. Per Pack Pcs
10001 33	1	0.01	0.010	360	12
10001 35	1	0.10	0.010	360	12
10001 41	2	0.02	0.015	360	12
10001 42	2	0.10	0.015	360	12
10001 53	5	0.05	0.045	360	12
10001 54	5	0.10	0.045	360	12
10001 60	10	0.10	0.075	360	12
10001 70	25	0.10	0.150	450	6



**Precision Graduated Pipettes, Clear Glass, Short-Line, Graduated to Tip, Amber Stain Graduation, (Serological type), Class 'B', EM-TECHCOLOR Brand** DIN EN ISO 835

Cat No.	Capacity ml	Graduation Division ml	Tolerance $\pm$ ml	Total Length $\pm$ 5mm	Qty. Per Pack Pcs
10001 29	0.5	0.010	0.0075	360	12



**Pipettes, Graduated, Clear Glass, Main Point Ring Graduation, Graduated to Tip, Amber Stain Graduation, (Serological type), Class 'AS', Conformity Certified, with Batch Certificate, Calibrated at 3 Graduation levels for accuracy and tolerance, EM-TECHCOLOR Brand** DIN EN ISO 835

Cat No.	Capacity ml	Graduation Division ml	Tolerance $\pm$ ml	Total Length $\pm$ 5mm	Qty. Per Pack Pcs
11001 33	1	0.01	0.007	360	12
11001 35	1	0.10	0.007	360	12
11001 41	2	0.02	0.010	360	12
11001 42	2	0.10	0.010	360	12
11001 53	5	0.05	0.030	360	12
11001 54	5	0.10	0.030	360	12
11001 60	10	0.10	0.050	360	12
11001 70	25	0.10	0.100	450	6



**Pipettes, Graduated, Clear Glass, Main Point Blue Ring Graduation, Graduated to Tip, Schellbach, (Serological type), Class 'AS', Conformity Certified, with Batch Certificate, Calibrated at 3 Graduation levels for accuracy and tolerance, EM-TECHCOLOR Brand** DIN EN ISO 835

Cat No.	Capacity ml	Graduation Division ml	Tolerance $\pm$ ml	Total Length $\pm$ 5mm	Qty. Per Pack Pcs
11101 33	1	0.01	0.007	360	12
11101 41	2	0.02	0.010	360	12
11101 53	5	0.05	0.030	360	12
11101 54	5	0.10	0.030	360	12
11101 60	10	0.10	0.050	360	12
11101 70	25	0.10	0.100	450	6



Pipettes, Graduated, Clear Glass, Short-Line, Graduated, Capacity Mark above shoulder, Amber Stain Graduation, (Mohr type), Class 'B', EM-TECHCOLOR Brand

DIN EN ISO 835

Cat No.	Capacity ml	Graduation Division ml	Tolerance $\pm$ ml	Total Length $\pm$ 5mm	Qty. Per Pack Pcs
10002 33	1	0.01	0.010	360	12
10002 35	1	0.10	0.010	360	12
10002 41	2	0.02	0.015	360	12
10002 42	2	0.10	0.015	360	12
10002 53	5	0.05	0.045	360	12
10002 54	5	0.10	0.045	360	12
10002 60	10	0.10	0.075	360	12
10002 70	25	0.10	0.150	450	6



Pipettes, Graduated, Clear Glass, Main Point Ring Graduation, Graduated Capacity Mark above Shoulder, Amber Stain Graduation, (Mohr type), Class 'AS', Conformity Certified, with Batch Certificate, Calibrated at 3 Graduation levels for accuracy and tolerance, EM-TECHCOLOR Brand

DIN EN ISO 835

Cat No.	Capacity ml	Graduation Division ml	Tolerance $\pm$ ml	Total Length $\pm$ 5mm	Qty. Per Pack Pcs
11002 33	1	0.01	0.007	360	12
11002 35	1	0.10	0.007	360	12
11002 41	2	0.02	0.010	360	12
11002 42	2	0.10	0.010	360	12
11002 53	5	0.05	0.030	360	12
11002 54	5	0.10	0.030	360	12
11002 60	10	0.10	0.050	360	12
11002 70	25	0.10	0.100	450	6



Pipette, Transfer, Volumetric, Accuracy as per Class 'A', Made from Heat Resistant, Low Expansion Borosilicate Glass, with Individual Certificate of Performance, RIVIERA

IS 1117: 1975

Cat No.	Capacity ml	Tolerance $\pm$ ml	Qty. Per Pack Pcs
72710 33	1	0.007	25
72710 40	2	0.010	25
72710 55	5	0.015	25
72710 60	10	0.020	25
72710 65	20	0.030	25
72710 70	25	0.030	25
72710 75	50	0.050	12
72710 80	100	0.080	6





Pipettes, Volumetric, Clear Glass, One Mark, **Blue** Inscription, Class 'AS',  
Conformity Certified with Batch Certificate, DWK

DIN EN ISO 835

Cat	Capacity	Tolerance	Total Length	Qty. Per Pack
No.	ml	± ml	± 10mm	Pcs
23339 0105	1	0.008	325	12
23339 0208	2	0.010	350	12
23339 0302	3	0.010	350	6
23339 0405	4	0.015	410	6
23339 0508	5	0.015	410	6
23339 0602	6	0.015	410	6
23339 0705	7	0.015	410	6
23339 0808	8	0.020	450	6
23339 0902	9	0.020	450	6
23339 1007	10	0.020	450	6
23339 1504	15	0.030	520	6
23339 2003	20	0.030	520	6
23339 2509	25	0.030	530	6
23339 5009	50	0.050	550	6



Pipettes, Volumetric, Clear Glass, One Mark, Amber Stain Graduation, Class 'AS',  
Conformity Certified with Batch Certificate, EM-TECHCOLOR Brand

DIN EN ISO 648

Cat	Capacity	Tolerance	Total Length	Qty. Per Pack
No.	ml	± ml	± 10mm	Pcs
13401 33	1	0.008	325	12
13401 40	2	0.010	350	12
13401 48	3	0.010	350	6
13401 50	4	0.015	410	6
13401 55	5	0.015	410	6
13401 56	6	0.015	410	6
13401 57	7	0.015	410	6
13401 58	8	0.020	450	6
13401 59	9	0.020	450	6
13401 60	10	0.020	450	6
13401 64	15	0.030	520	6
13401 65	20	0.030	520	6
13401 70	25	0.030	530	6
13401 75	50	0.050	550	6



Pipettes, Volumetric, Clear Glass, Two Marks, Amber Stain Graduation, Class 'AS',  
Conformity Certified with Batch Certificate, EM-TECHCOLOR Brand

DIN EN ISO 648

Cat	Capacity	Tolerance	Total Length	Qty. Per Pack
No.	ml	± ml	± 10mm	Pcs
13402 33	1	0.008	325	12
13402 40	2	0.010	350	12
13402 48	3	0.010	350	6
13402 55	5	0.015	410	6
13402 60	10	0.020	450	6
13402 64	15	0.030	520	6
13402 65	20	0.030	520	6
13402 70	25	0.030	530	6
13402 75	50	0.050	550	6



Pipettes, Transfer, Volumetric, Accuracy as per Class 'B',  
Made from Heat Resistant, Low Expansion Borosilicate Glass, RIVIERA

IS 1117 : 1975



Cat	Capacity	Tolerance	Qty. Per Pack
No.	ml	± ml	Pcs
72700 33	1	0.015	25
72700 40	2	0.020	25
72700 55	5	0.030	25
72700 60	10	0.040	25
72700 65	20	0.060	25
72700 70	25	0.060	25
72700 75	50	0.100	12
72700 80	100	0.160	6

Culture Tubes, Round Bottom, with DIN Thread with Screw-Cap made of PP and Rubber Gasket, (Suitable for food use) Sterilizable at 1 bar 121°C, AR-GLAS



Cat	O. D. x Length	DIN Thread	Wall Thickness	Qty. Per Pack
No.	mm	GL	mm	Pcs
23175 115	12 x 100	14	1	100
23175 145	16 x 100	18	1	100
23175 215	16 x 160	18	1	100
23175 235	18 x 180	18	1	100

Supplies will be without any brand mark on the Culture Tubes.

Culture Tubes, Round Bottom, with Screw Cap and Silicone Liner, Made from Heat Resistant, Low Expansion Borosilicate Glass, RIVIERA



Cat	Capacity	O. D. x Length	Qty. Per Pack
No.	ml	mm	Pcs
71300 055	5	16 x 75	100
71300 105	10	16 x 125	100
71300 305	30	25 x 100	100
71300 605	60	25 x 200	100

Culture Tubes, Round Bottom, AMBER, with Screw Cap and Silicone Liner, Made from Heat Resistant, Low Expansion Borosilicate Glass, RIVIERA



Cat	Capacity	O. D. x Length	Qty. Per Pack
No.	ml	mm	Pcs
71310 055	5	16 x 75	100
71310 105	10	16 x 125	100
71310 305	30	25 x 100	100
71310 605	60	25 x 200	100