



### Calibration Solutions \*

Wide range of calibration solutions for pH, ORP, Conductivity, TDS and Salinity. Calibration solutions come in a convenient ready to use 480 ml leak-proof bottles. Calibration must be done on regular basis for maintaining accuracy of the instrument. It also helps negate the effect of electrode ageing and contamination. Minimum shelf life of the solutions is as mentioned.

#### pH Buffer Solutions ( 480 ml)

Cat. No.	Product Code No.	Description	Shelf Life from Date of Packing
8005 001	ECBU4BT	pH 4.01 buffer	36 months
8005 061	ECBU4BTC	pH 4.01 buffer (colored)	36 months
8005 002	ECBU7BT	pH 7.00 buffer	36 months
8005 062	ECBU7BTC	pH 7.00 buffer (colored)	36 months
8005 003	ECBU9BT	pH 9.00 buffer	36 months
8005 004	ECBU10BT	pH 10.01 buffer	36 months
8005 064	ECBU10BTC	pH 10.01 buffer (colored)	36 months
8005 005	ECBU12BT	pH 12.45 buffer	36 months
8005 006	ECBU1BT	pH 1.68 buffer	36 months
8005 007	ECBU686BT	pH 6.86 buffer	36 months
8005 008	ECBU918BT	pH 9.18 buffer	36 months

#### Conductivity Solutions ( 480 ml)

Cat. No.	Product Code No.	Description	Shelf Life from Date of Packing
8006 008	ECCON84BT	84 $\mu$ S KCL	12 months
8006 001	ECCON100BT	100 $\mu$ S KCL	12 months
8006 002	ECCON500BT	500 $\mu$ S KCL	36 months
8006 003	ECCON1413BT	1413 $\mu$ S KCL	36 months
8006 004	ECCON2764BT	2764 $\mu$ S KCL	36 months
8006 005	ECCON1288BT	12.88 mS KCL	36 months
8006 006	ECCON1118BT	111.8 mS KCL	36 months
8006 007	ECCON5000BT	5.0 mS KCL	36 months

#### TDS Solutions ( 480 ml)

Cat. No.	Product Code No.	Description	Shelf Life from Date of Packing
8007 004	EC44250BT	50 ppm	12 months
8007 001	EC442300BT	300 ppm	36 months
8007 002	EC4421000BT	1000 ppm	36 months
8007 003	EC4423000BT	3000 ppm	36 months

#### ORP Calibration Standards

Cat. No.	Product Code No.	Description	Shelf Life from Date of Packing
8047 001	ECORPPRE	Pretreatment (475 mV)	36 months
8047 002	ECORPQUIN	Quinhydrone 263 (263 mV, $\pm$ 30 mV, 25 °C)	36 months
8047 003	ECORPQUIN086	Quinhydrone 86 (86 mV)	36 months

#### Salinity Calibration Standards

Cat. No.	Product Code No.	Description	Shelf Life from Date of Packing
8048 001	ECNACL5PPT	5 ppt NaCl	36 months
8048 002	ECNACL25PPT	25 ppt NaCl	36 months
8048 003	ECNACL45PPT	45 ppt NaCl	36 months