

Orbital Shaker - Heavy Duty, Rivotek



Specifications

- Body Material - MS Powder coated.
- Tray Material – Stainless steel SS 304
- Shaking Speed: 20 – 250 Orbital strokes per min (Without Load)
- Operates on 230 Volts AC, Single Phase

Special Features

- **VFD Drive:** (Variable Frequency Drive) of 230 V, I Phase Input.
- **AC motor:** Brushless Induction Motor, for 24 hrs operation.
- **Digital Voltmeter:** With LCD display for better view in the dark room.

Cat No.	Tray size in mm	Digital Speed Indication	Shaking Speed (without load) in RPM
50081 0750	750 x 750	Yes	20 - 250



Ultrasonic Cleaners, Rivotek

Specifications :

- Bath Material : SS 304
- Outer Body Material & Lid : SS 304
- Ultrasonic Generator : Advanced latest MOSFET/IGBT based SMPS
- Frequency : 30 ± 3 KHz
- Transducer : Imported make PZT crystals bonded on the base of the tank with weld bond technique.
- SS Basket : Provided
- Built-in Digital Timer : Provided
- Electric supply : 230V/ 50Hz AC, single phase.
- Drain : Drain valve provided at the rear side of the tank (not included for 1.5 litre model)

Cat No.	Tank Cap in litres	Tank Clear Size in mm	Ultrasonic Power Watts
1900 0015	1.5	240 x 135 x 65	50
1900 0055	5.0	230 x 230 x 100	100
1900 0090	9.0	300 x 240 x 150	200
1900 0020	20.0	350 x 300 x 200	400
1900 0030	30.0	375 x 375 x 200	500
1900 0050	50.0	550 x 300 x 300	750