




## GEL ROCKER



The Rivotek Gel Rocker can be used for many applications like gel staining, culture aeration, bathing, blotting membranes etc. A metal platform with non slip rubber mat moves in a gentle 'See Saw' motion. The movement is controlled precisely with a DC motor. Cushion and rubber feet prevent unit from 'slipping'. Time can be set from 1min to 99 mins. or continuous.

### Specification

<b>Model :</b>	<b>50051 003</b>
Motor :	D.C.
Maximum RPM :	60
Speed controller :	Analogue
Timer :	Digital
Time range :	1-99mins/Cont.
Platform Size :	25 x 25 cm
Weight :	4 kgs.
Dimension :	25 x 25 x 17.5 cms
Input Voltage :	230V, 50Hz, A.C. Supply

-----  
**Riviera Glass Pvt. Ltd.** ♦C-2116, Oberoi Garden Estates, Chandivali, Andheri (E), Mumbai – 400072 ♦ □Tel.No. 022-28475228 ♦ Email id- info@rivieraglass.co.in.,  
rivieraglass1995@gmail.com | [www.rivieraglass.co.in](http://www.rivieraglass.co.in) ♦