



B.O.D. INCUBATOR

Specifications:

- Material of Construction: Interior cabinet- SS 304 Grade
- Exterior cabinet: MS Powder coated.
- Temperature range & Accuracy : 5 °C to 60 °C & ± 0.2 °C
- Temperature Sensor: Class "A", PT 100 RTD
- Display/ Control: Digital Microprocessor based Temperature Controller with resolution 0.1 °C and accuracy ± 0.2 °C with non-volatile storage memory for the set Point.
- Refrigeration: CFC free compressor utilizing R134a Eco-friendly refrigeration
- Observation: Inner transparent acrylic door and outer door key lockable.
- Illumination: Florescent light with door switch
- Feet: Castor wheels
- Mains Power: 230V/ 50hz AC



Ordering Information:

Cat. No	Capacity Ltrs./Cu Ft.	Internal Dimensions W x D x H (mm)	No. of Trays
6101 0004	120 / 04	45 x 45 x 60	2
6101 0008	200 / 08	60 x 60 x 60	3
6101 0012	324 / 12	60 x 60 x 90	3
6101 0016	440 / 16	60 x 60 x 125	4
6101 0021	600/21	80 x 60 x 125	5
6101 0035	1000 / 35	100 x 80 x 125	5