




---

**Optional Accesories for Extraction and Kjeldahl Apparatus, Rivotek**


---

Cat	Description
No.	
9002 0001	Finger Clamp with Boss Head

---


**Test Tube Rotator, Rivotek**

Speed Variable from 5 - 50 RPM (Without Load)  
 Permanent magnet DC geared motor.  
 Noise-free silent operation.  
 Constant current speed control circuit for steady speed.  
 Powder coated M.S.body.

**(Note : While placing the order, please specify the dimensions of the tubes you propose to use as per details given below)**

---

Cat	Description
No.	
50022 001	Test Tube Rotator with 0.5 ml Eppendorf tubes 24 Nos. & 1.5 ml Eppendorf tubes 24 Nos.
50022 002	Test Tube Rotator with $\varnothing$ 16/19 mm Test tubes 12 Nos.
50022 003	Test Tube Rotator with $\varnothing$ 30 mm Test tubes 8 Nos.

---

**Optional Accesories for Test Tube Rotator, RIVOTEK**


---

Cat	Description
No.	
50022 501	Digital Timer (0 - 99 min)
50022 502	Digital RPM

---


**Universal Cell-Mixer, Plate Type, Rivotek**

Suitable for mini tissue culture work and blood suspensions without sedimentation or foaming.  
 Also suitable for vigorous agitation and end to end mixing.  
 Tubes can be mounted at any angle from horizontal to vertical, tubes with respect to axis of rotation.  
 Different sizes from 12 mm to 25 mm dia can be accommodated.  
 Speed Variable 5 - 50 rpm (wihout load)

---

Cat	Description
No.	
50023 001	Universal Cell Mixer with 0.5 ml ( $\varnothing$ 8 mm) Eppendorf tubes 12 Nos.
50023 002	Universal Cell Mixer with 1.5 ml ( $\varnothing$ 11 mm) Eppendorf tubes 12 Nos.
50023 003	Universal Cell Mixer with $\varnothing$ 13 to 27 mm Test tubes 8 Nos.

---