

Compulsory Accessories

Spring loaded SS lotus clamps (removable) for holding the Flasks in the Water Shakers to be ordered separately

Cat No.	Flask capacity in ml	No. of SS Lotus Clamp required for tray size		
		250 x 250 mm	325 x 325 mm	350 x 350 mm
50081 501	50	9	16	20
50081 502	100	9	16	20
50081 503	250	5	8	9
50081 504	500	4	8	6
50081 505	1000	-	4	-

Optional Accesories for Water Bath, RIVOTEK

Cat No.	Description
50022 503	Digital Safety Controller
50022 502	Digital RPM


Microplate Shaker, Rivotek

- Body Material is MS painted with Epoxy Powder. Tray Material is SS.
- Shaking speed: 100-1500 strokes per min
- Permanent Magnet Geared DC motor with 1mm dia. Orbital motion is used.

Cat No.	Tray size
50091 001	180 x 150 mm. With springs 2 multi well plates can be placed on the tray

Water Bath Double-walled (Thermostatic Controlled)


Specifications

- Inner chamber is made of thick gauge S.S. 304 grade sheet, with minimum joint. STD/GMP
- Outer MS with Powder Coating only for - STD model
- Outer made up of S.S. 304 grade sheet only for GMP model
- Specially made extra long heating elements for high heat dissipation, uniformity & longer life.
- The gap between the two walls is filled with glass wool.
- Lid is made up of S.S. 304 with 6 holes & 12 holes of 3" dia.
- Provided with indicating lamps and on/off switch.
- Temp range is above (Ambient to +5 to 99 °C)
- Operation on 240 volts AC

Cat No.	Size in inches	Holes	Watt	Inner Chamber	Outer Chamber
50192 000	14 x 10 x 4	6	1000	S. S. 304	MS Powder Coated
50192 001	17 x 13 x 4	12	2000	S. S. 304	MS Powder Coated
50191 000	14 x 10 x 4	6	1000	S. S. 304	S. S. 304
50191 001	17 x 13 x 4	12	2000	S. S. 304	S. S. 304
