




**Bottles, Culture, Roux-Type, Conical Neck, Eccentric, DURAN**



Cat No.	Capacity ml	Breadth mm	Depth mm	Height mm	Qty. Per Pack Pcs
21571 58	1200	132	56	275	10

**Bottles, Wash, LDPE Plastic, RIVIERA**

Bottles & delivery tubes made of chemically resistant low density polyethylene, flexible with screw cap.



Cat No.	Capacity ml	Qty. Per Pack Pcs
74200 36	250	10
74200 44	500	10

**Burettes, Straight Bore Stopcock, Accuracy as per Class 'A', with Individual Certificate of Performance, Made from Heat Resistant, Low Expansion Borosilicate Glass, RIVIERA** IS 1997 : 2008



Cat No.	Capacity ml	Graduation Division ml	Tolerance ±ml	App. Total Height mm	Qty. Per Pack Pcs
70301 60	10	0.05	0.02	530	5
70301 70	25	0.10	0.05	530	5
70301 75	50	0.10	0.05	780	5
70301 80	100	0.10	0.10	830	5

**Burettes, Rotaflow Needle Valve Stopcock, Accuracy as per Class 'A', with Individual Certificate of Performance, Made from Heat Resistant, Low Expansion Borosilicate Glass, RIVIERA** IS 1997 : 2008



Cat No.	Capacity ml	Graduation Division ml	Tolerance ±ml	App. Total Height mm	Qty. Per Pack Pcs
70401 60	10	0.05	0.02	530	5
70401 70	25	0.10	0.05	530	5
70401 75	50	0.10	0.05	780	5
70401 80	100	0.10	0.10	830	5

**Burettes, Schellbach, Main Point Ring Graduation, Straight Stopcock with Glasskey, Blue Graduation, Class 'AS', Conformity Certified, with Batch Certificate, Calibrated at 3 Graduation levels for accuracy and tolerance, EM-TECHCOLOR Brand** DIN EN ISO 385



Cat No.	Capacity ml	Graduation Division ml	Tolerance ±ml	Qty. Per Pack Pcs
31401 60	10	0.02	0.02	2
31401 70	25	0.05	0.03	2
31401 75	50	0.10	0.05	2



**Burettes Straight Bore Stopcock, Accuracy as per Class 'B',** IS 1997 : 2008  
**Made from Heat Resistant, Low Expansion Borosilicate Glass, RIVIERA**

Cat	Capacity	Graduation	Tolerance	App. Total	Qty. Per
No.	ml	Division	$\pm$ ml	Height	Pack
		ml		mm	Pcs
70300 60	10	0.05	0.05	530	5
70300 70	25	0.10	0.10	530	5
70300 75	50	0.10	0.10	780	5
70300 80	100	0.10	0.20	830	5



**Burettes, Rotaflow Needle Valve Stopcock, Accuracy as per Class 'B',** IS 1997 : 2008  
**Made from Heat Resistant, Low Expansion Borosilicate Glass, RIVIERA**

Cat	Capacity	Graduation	Tolerance	App. Total	Qty. Per
No.	ml	Division	$\pm$ ml	Height	Pack
		ml		mm	Pcs
70400 60	10	0.05	0.05	530	5
70400 70	25	0.10	0.10	530	5
70400 75	50	0.10	0.10	780	5
70400 80	100	0.10	0.20	830	5



**Burettes, Schellbach, Main Point Ring Graduation, Straight Valve Stopcock with PTFE Spindle,** DIN EN ISO 385  
**Blue Graduation, Class 'AS', Conformity Certified, with Batch Certificate,**  
**Calibrated at 3 Graduation levels for accuracy and tolerance, EM-TECHCOLOR Brand**

Cat	Capacity	Graduation	Tolerance	Qty. Per
No.	ml	Division	$\pm$ ml	Pack
		ml		Pcs
31402 60	10	0.02	0.02	2
31402 70	25	0.05	0.03	2
31402 75	50	0.10	0.05	2



**Burettes, Automatic, Straight Bore Stopcock with Glasskey, White Graduation,** IS 1997 : 2008  
**Accuracy as per Class 'A', with Individual Certificate of Performance, Complete with Reservoir,**  
**Made from Heat Resistant, Low Expansion Borosilicate Glass, RIVIERA**

Cat	Capacity	Qty. Per
No.	ml	Pack
		Pcs
70312 60	10	1
70312 70	25	1
70312 75	50	1
70312 80	100	1



**Burettes, Automatic, Rotaflow Stopcock, White Graduation, Accuracy as per Class 'A',** IS 1997 : 2008  
**with Individual Certificate of Performance, Complete with Reservoir,**  
**Made from Heat Resistant, Low Expansion Borosilicate Glass, RIVIERA**

Cat	Capacity	Qty. Per
No.	ml	Pack
		Pcs
70412 60	10	1
70412 70	25	1
70412 75	50	1
70412 80	100	1



**Burettes, Automatic, Mounted on Reservoir, Schellbach, Main Point Ring Graduation, Straight Stopcock with Glasskey, Blue Graduation, Class 'AS', Conformity Certified with Batch Certificate, Calibrated at 3 Graduation levels for accuracy and tolerance, Made from EM-TECHCOLOR Brand**  
Batch Certificate shall be for Burette only (Cat. No. 31401..)

Cat	Capacity	Qty. Per Pack
No.	ml	Pcs
314001 60	10	1
314001 70	25	1
314001 75	50	1

\* Reservoir Capacity : 500 ml for 10 ml and 2000 ml for 25 & 50 ml size.



**Burette, Automatic, Mounted on Reservoir, Schellbach, Main Point Ring Graduation, Straight Valve Stopcock with PTFE Spindle, Blue Graduation, Class 'AS', Conformity Certified with Batch Certificate, Calibrated at 3 Graduation levels for accuracy and tolerance, Made from EM-TECHCOLOR Brand**  
Batch Certificate shall be for Burette only (Cat. No. 31402..)

Cat	Capacity	Qty. Per Pack
No.	ml	Pcs
314002 60	10	1
314002 70	25	1
314002 75	50	1

\* Reservoir Capacity : 500 ml for 10 ml and 2000 ml for 25 & 50 ml size.



**Burettes, Automatic, Straight Bore Stopcock with Glasskey, White Graduation, Accuracy as per Class 'B', Complete with Reservoir, Made from Heat Resistant, Low Expansion Borosilicate Glass, RIVIERA** IS 1997: 2008

Cat	Capacity	Qty. Per Pack
No.	ml	Pcs
70311 60	10	1
70311 70	25	1
70311 75	50	1

\* Reservoir Capacity : 500 ml for 10 ml and 2000 ml for 25 & 50 ml size.



**Burettes, Automatic Zero, Rotaflow Stopcock, White Graduation, Accuracy as per Class 'B', Complete with Reservoir, Made from Heat Resistant, Low Expansion Borosilicate Glass, RIVIERA** IS 1997: 2008

Cat	Capacity	Qty. Per Pack
No.	ml	Pcs
70411 60	10	1
70411 70	25	1
70411 75	50	1

\* Reservoir Capacity : 500 ml for 10 ml and 2000 ml for 25 & 50 ml size.



**Reservoir for Automatic Burette, Made from DURAN**

Cat	Capacity	Joint Size	Qty. Per Pack
No.	ml		Pcs
2180 144	500	B 29	1
2180 163	2000	B 34	1
2180 173	5000	B 34	1

