



**Beakers, Low Form, with Graduation and Spout,
Made from Heat Resistant, Low Expansion Borosilicate Glass, RIVIERA**



Cat	Capacity	O. D.	Height	Qty. Per Pack
No.	ml	mm	mm	Pcs
70100 08	5	22	30	10
70100 10	10	26	35	10
70100 14	25	34	50	10
70100 17	50	42	60	10
70100 24	100	50	70	10
70100 29	150	60	80	10
70100 36	250	70	95	12
70100 41	400	80	110	12
70100 44	500	83	115	6
70100 48	600	90	125	6
70100 54	1000	105	145	6
70100 63	2000	130	185	4
70100 68	3000	150	210	1
70100 73	5000	170	270	1
70100 86	10000	220	350	1

*5 ml & 10 ml will be supplied without graduation

Beakers, Low Form, with Graduation and Spout, DURAN **DIN 12331, ISO 3819, USP Standard**



Cat	Capacity	O. D.	Height	Qty. Per Pack
No.	ml	mm	mm	Pcs
21106 14	25	34	50	10
21106 17	50	42	60	10
21106 24	100	50	70	10
21106 29	150	60	80	10
21106 36	250	70	95	10
21106 41	400	80	110	10
21106 48	500/600	90	125	10
21106 54	1000	105	145	10
21106 63	2000	132	185	10
21106 68	3000	152	210	4
21106 73	5000	170	270	1
21106 86	10000	217	350	1

*10000 ml Capacity is a non - DIN non - ISO size and will be supplied without graduation

**Beakers, Tall Form, with Graduation and Spout,
Made from Heat Resistant, Low Expansion Borosilicate Glass, RIVIERA**



Cat	Capacity	O. D.	Height	Qty. Per Pack
No.	ml	mm	mm	Pcs
70110 17	50	38	70	24
70110 24	100	48	80	24
70110 29	150	54	95	12
70110 36	250	60	120	12
70110 44	500	75	140	6
70110 54	1000	95	180	6
70110 63	2000	120	240	1



Beakers, Tall Form, with Graduation and Spout, DURAN				DIN 12331, ISO 3819, USP Standard
Cat	Capacity	O. D.	Height	Qty. Per Pack
No.	ml	mm	mm	Pcs
21116 17	50	38	70	10
21116 24	100	48	80	10
21116 29	150	54	95	10
21116 36	250	60	120	10
21116 41	400	70	130	10
21116 48	500/600	80	150	10
21116 54	1000	95	180	10
21116 63	2000	120	240	10
21116 68	3000	135	280	2



Beakers, Tall Form, with Graduation without Spout, Made from Heat Resistant, Low Expansion Borosilicate Glass, RIVIERA				
Cat	Capacity	O. D.	Height	Qty. Per Pack
No.	ml	mm	mm	Pcs
70120 36	250	60	120	12
70120 44	500	75	140	6
70120 54	1000	95	180	6



Beakers, Tall Form, with Graduation without Spout, DURAN				USP Standard
Cat	Capacity	O. D.	Height	Qty. Per Pack
No.	ml	mm	mm	Pcs
21117 17	50	38	70	10
21117 24	100	48	80	10
21117 29	150	54	95	10
21117 36	250	60	120	10
21117 41	400	70	130	10
21117 48	500/600	80	150	10
21117 54	1000	95	180	10



Beakers, Philips with Spout, Made from Heat Resistant, Low Expansion Borosilicate Glass, RIVIERA				
Cat	Capacity	Base Dia	Height	Qty. Per Pack
No.	ml	mm	mm	Pcs
70130 36	250	75	113	10
70130 44	500	90	170	10



Beakers, Philips with Spout, DURAN					
Cat	Capacity	Base Dia	Mouth Dia	Height	Qty. Per Pack
No.	ml	mm	mm	mm	Pcs
21141 29	150	59	49	87	10
21141 36	250	68	52	105	10
21141 44	500	86	59	142	10