

**OPTIONAL ACCESSORIES FOR DISTILLATION UNITS**

Cat No.	Description	Qty. Per Pack Pcs
72250 000	Auto Cut-Off System with Audio Alarm for single stage (Capacities 1 to 5 ltrs/hr) and also for double stage (Capacities upto 2.5 ltrs/hr only)	1
72250 100	Auto Cut-Off System with Audio Alarm for double stage (Capacities 3 ltrs/hr and above only)	1
72260 000	Water Softner for Boiler Feed	1
72260 001	Resin for Water Softner	1
72260 500	DM Plant	1

**Warranty Clause for Riviera Distillation Units:**

The above equipment is warranted with following conditions:

- The items supplied will carry a warranty of 12 months from the date of the invoice against manufacturing defects. The company will repair or replace the parts going defective free of cost during the warranty period. Vacuumatic items and consumables like bulb, fuse, heater, battery, etc. are not covered under warranty.
- Glass parts do not carry any warranty after installation.
- Complaints arising out of mishandling willful or through ignorance, negligence or resulting from accidents, attempted repairs, and inadequate power supply, fluctuation in voltage, general wear and tear are not covered under this warranty.
- Complaints should be communicated to the warrantor in writing.

**All Glass Filter Flask for Filter Holder 47mm**

Excellent suitability for vacuum filtration of liquids for analysis, especially suitable for HPLC set up, both solvent and mobile phases. Very good compatability with chemicals including very aggressive solvents - all parts made of borosilicate glass.

**Specifications :**

Filter size : 47mm dia.

Effective Filtration area : 9.60 cm<sup>2</sup>

Filter flask capacity : 1000 ml

Funnel capacity : 300 ml

Spring clamp : Anodised Aluminium

Pressure : Vacuum only

Cat No.	Description	Qty. Per Pack Pcs
73200 544	Complete set	1
	<b>Parts only</b>	
73200 504	Flask for all glass filter holder 47mm	1
73200 502	Funnel for all glass filter holder 47mm	1
73200 503	Filter base for all glass filter holder 47mm	1
73200 501	Clamp	1

