



**Flasks, Erlenmeyer, Conical, Narrow Mouth, with Graduation
Made from Heat Resistant, Low Expansion Borosilicate Glass, RIVIERA**

Cat No.	Capacity ml	Dia Base mm	Neck Dia mm	Height mm	Qty. Per Pack Pcs
71400 10	10	37	15	60	10
71400 14	25	42	22	75	10
71400 17	50	51	22	90	10
71400 24	100	64	22	105	10
71400 29	150	72	28	124	10
71400 36	250	85	34	145	10
71400 44	500	105	34	180	10
71400 54	1000	131	42	220	10
71400 63	2000	166	50	280	4
71400 68	3000	187	50	310	1
71400 73	5000	220	50	365	1



Flasks, Erlenmeyer, Conical, Narrow Mouth, with Graduation, DURAN DIN 12380, ISO 1773

Cat No.	Capacity ml	Dia Base mm	Neck Dia mm	Height mm	Qty. Per Pack Pcs
21216 14	25	42	22	75	10
21216 17	50	51	22	90	10
21216 24	100	64	22	105	10
21990 27	150	74	28	118	10
21216 36	250	85	34	145	10
21216 44	500	105	34	180	10
21216 54	1000	131	42	220	10
21216 63	2000	166	50	280	10
21216 68	3000	187	51	310	2
21216 73	5000	220	51	365	1




Flasks, Erlenmeyer, Conical with PP Screw Cap, Autoclavable up to 140°C, Made from DURAN

Cat No.	Capacity ml	Qty. Per Pack Pcs
21803 365	250	10
21803 445	500	10
21803 545	1000	10



**Flasks, Iodine Determination with Glass Stopper,
Made from Heat Resistant, Low Expansion Borosilicate Glass, RIVIERA**

Cat No.	Capacity ml	Joint Size	Qty. Per Pack Pcs
74400 24	100	B 24	10
74400 36	250	B 24	10
74400 44	500	B 24	10
74400 54	1000	B 24	10



**Flasks, Erlenmeyer, Conical, Wide Mouth,
Made from Heat Resistant, Low Expansion Borosilicate Glass, RIVIERA**

Cat No.	Capacity ml	Dia Base mm	Neck Dia mm	Height mm	Qty. Per Pack Pcs
71410 24	100	70	34	110	24
71410 36	250	85	50	142	6
71410 44	500	105	50	175	6
71410 54	1000	131	56	220	6



**Flasks, Erlenmeyer, Conical, Wide Mouth,
with Graduation, DURAN** DIN 12385

Cat No.	Capacity ml	Dia Base mm	Neck Dia mm	Height mm	Qty. Per Pack Pcs
21226 14	25	43	31	70	10
21226 17	50	51	34	85	10
21226 24	100	64	34	105	10
21226 36	250	85	50	140	10
21226 44	500	105	50	175	10
21226 54	1000	131	50	220	10
21226 63	2000	153	72	276	10




**Flasks, Flat Bottom, Narrow Neck,
Made from Heat Resistant, Low Expansion Borosilicate Glass, RIVIERA**

Cat No.	Capacity ml	Max. Dia mm	Neck Dia mm	Height mm	Qty. Per Pack Pcs
73900 14	25	42	20	90	10
73900 17	50	51	26	90	10
73900 24	100	64	26	105	10
73900 29	150	75	28	132	10
73900 36	250	85	34	140	10
73900 44	500	105	34	165	10
73900 54	1000	131	42	190	10
73900 64	2000	166	50	230	6
73900 68	3000	186	50	300	1
73900 73	5000	224	58	320	1
73900 86	10000	279	65	360	1



Flasks, Flat Bottom, Narrow Neck, with Beaded Rim, DURAN DIN 12347, ISO 1773

Cat No.	Capacity ml	Max. Dia mm	Neck Dia mm	Height mm	Qty. Per Pack Pcs
21711 36	250	85	35	140	10
21711 44	500	105	34	171	10
21711 54	1000	131	42	200	10
21711 64	2000	166	42	250	10
21711 68	3000	185	50	250	1
21711 73	5000	223	50	290	1



Flasks, Pear Shape, Suitable for Rotary Evaporators with Ground Socket 29/32, DURAN

Cat No.	Capacity ml	Max Dia mm	Height mm	Qty. Per Pack Pcs
24120 37	250	81	140	10
24120 46	500	101	170	10
24120 56	1000	126	210	10



**Flasks, Round Bottom, Narrow Neck,
Made from Heat Resistant, Low Expansion Borosilicate Glass, RIVIERA**

Cat No.	Capacity ml	Max. Dia mm	Neck Dia mm	Height mm	Qty. Per Pack Pcs
74000 14	25	42	20	90	10
74000 17	50	51	26	95	10
74000 24	100	64	26	110	10
74000 29	150	75	28	137	10
74000 36	250	85	34	143	10
74000 44	500	105	34	168	10
74000 54	1000	131	42	200	10
74000 64	2000	166	51	240	6
74000 68	3000	190	50	305	1
74000 73	5000	224	58	350	1
74000 86	10000	279	65	380	1
74000 91	20000	350	78	505	1



**Flasks, Distilling with Side Arm,
Made from Heat Resistant, Low Expansion Borosilicate Glass, RIVIERA
Flasks of capacity 250 ml onwards made from DURAN**

Cat No.	Capacity ml	Approximate O. D. x Height mm	Approximate Neck O. D. mm	Qty. Per Pack Pcs
22253 24	100	64 x 210	22	5
22253 25	125	68 x 214	19	5
22253 36	250	88 x 226	24	2
22253 44	500	110 x 250	28	2
22253 54	1000	130 x 305	30	1
22253 64	2000	170 x 318	40	1

Flasks, Kjeldahl, Made from Heat Resistant, Low Expansion Borosilicate Glass, RIVIERA

Cat No.	Capacity ml	Max. Dia mm	Neck Dia mm	Height mm	Qty. Per Pack Pcs
74500 44	500	100	29	340	6
74500 51	800	118	34	365	6



Flasks, Kjeldahl, with I/C Joint, Made from Heat Resistant, Low Expansion Borosilicate Glass, RIVIERA

Cat No.	Capacity ml	Joint Size	Qty. Per Pack Pcs
74501 44	500	B 24	6
74501 51	800	B 24	6



Flasks, Kjeldahl, DURAN

DIN 12360

Cat No.	Capacity ml	Max. Dia mm	Neck Dia mm	Height mm	Qty. Per Pack Pcs
21231 44	500	101	34	300	10
21231 51	750	115	34	340	10



Flasks, Kjeldahl, with I/C Joint, Made from DURAN

Cat No.	Capacity ml	Joint Size	Qty. Per Pack Pcs
217030 44	500	B 24	10
217030 51	750	B 24	10

