

<b>Stopcock with PTFE Spindle for Vacuum Desiccator</b>			
<b>Cat</b>	<b>Joint Size</b>	<b>Length</b>	<b>Qty. Per Pack</b>
<b>No.</b>		<b>mm</b>	<b>Pcs</b>
24799 04	B 24	124	1



<b>Dishes, Crystallizing, Made from Heat Resistant, Low Expansion Borosilicate Glass, RIVIERA</b>			
<b>Cat</b>	<b>Capacity</b>	<b>O. D. x Height</b>	<b>Qty. Per Pack</b>
<b>No.</b>	<b>ml</b>	<b>mm</b>	<b>Pcs</b>
71600 44	300	95 x 55	6
71600 54	900	140 x 75	4
71600 59	2000	190 x 90	1



<b>Dishes, Crystallizing, with Spout, DURAN</b>			<b>DIN 12338</b>
<b>Cat</b>	<b>Capacity</b>	<b>O. D. x Height</b>	<b>Qty. Per Pack</b>
<b>No.</b>	<b>ml</b>	<b>mm</b>	<b>Pcs</b>
21311 44	300	95 x 55	10
21311 54	900	140 x 75	10
21311 59	2000	190 x 90	10



<b>Dishes, Evaporating, Made from Heat Resistant, Low Expansion Borosilicate Glass, RIVIERA</b>			
<b>Cat</b>	<b>Capacity</b>	<b>Max. O. D. x Height</b>	<b>Qty. Per Pack</b>
<b>No.</b>	<b>ml</b>	<b>mm</b>	<b>Pcs</b>
71700 41	90	80 x 45	10
71700 44	170	105 x 55	2
71700 54	600	150 x 80	2



<b>Dishes, Evaporating, with Spout, DURAN</b>			<b>DIN 12336</b>
<b>Cat</b>	<b>Capacity</b>	<b>Max. O. D. x Height</b>	<b>Qty. Per Pack</b>
<b>No.</b>	<b>ml</b>	<b>mm</b>	<b>Pcs</b>
21301 41	90	80 x 45	10
21301 44	170	95 x 55	10
21301 54	600	140 x 80	10



<b>Dishes, Petri, Made from Heat Resistant, Low Expansion Borosilicate Glass, RIVIERA</b>			
<b>Cat</b>	<b>O. D. x Height</b>		<b>Qty. Per Pack</b>
<b>No.</b>	<b>mm</b>		<b>Pcs</b>
71500 40	50 x 17		10
71500 42	80 x 17		10
71500 48	100 x 17		10
71500 52	150 x 20		10
71500 59	200 x 20		10

