

Riviera

REAGENT BOTTLE

Narrow Mouth with PP Screw Cap



REAGENT BOTTLES WITH PP SCREW CAPS, RIVIERA

The Ideal bottles for storage and sampling

RIVIERA Reagent Bottles, made from 3.3 Borosilicate Glass, narrow neck, with PP Screw Caps are now available in both Clear and Amber stained glass. Highly suitable for pouring of liquids as the necks have a PP pouring ring which prevents the liquids from dripping. Bottles with PP Screw Caps and Pouring Rings are fully autoclavable up to 140°C

Amber stained bottles are particularly suitable for storage of light sensitive chemicals

Both the bottles, be it clear or amber have the features

normally associated with borosilicate glass 3.3 expansion - chemically neutral, heat resistant and thermal shock proof

Bottles from 100ml to 20,000ml have single thread size with GL 45 caps. For 25ml and 50ml the thread and cap sizes are GL 25 and GL 32 respectively

25ml Bottles, Being small in size does not come with Pouring Rings

While autoclaving do remember to close the bottles loosely to ensure pressure equalisation.



REAGENT BOTTLE, Narrow Mouth with PP Screw Cap, Clear, RIVIERA

Cat No.	Capacity ml	Thread GL	O. D mm	Height mm	Qty Per Pack
71801 145	25	25	36	70	10
71801 175	50	32	46	87	10
71801 245	100	45	56	100	10
71801 365	250	45	70	138	10
71801 445	500	45	86	177	10
71801 545	1000	45	101	225	10
71801 635	2000	45	136	260	6
71801 695	3000	45	158	295	1
71801 735	5000	45	181	320	1
71801 865	10000	45	235	405	1
71801 915	20000	45	296	505	1

25 ml Bottles will be supplied without pouring ring.
Height indicated is of Bottle Only



REAGENT BOTTLE, Narrow Mouth with PP Screw Cap, Amber, RIVIERA

Cat No.	Capacity ml	Thread GL	O. D mm	Height mm	Qty Per Pack
71806 145	25	25	36	70	10
71806 175	50	32	46	87	10
71806 245	100	45	56	100	10
71806 365	250	45	70	138	10
71806 445	500	45	86	177	10
71806 545	1000	45	101	225	10
71806 635	2000	45	136	260	6
71806 695	3000	45	158	295	1
71806 735	5000	45	181	320	1
71806 865	10000	45	235	405	1
71806 915	20000	45	296	505	1

25 ml Bottles will be supplied without pouring ring.
Height indicated is of Bottle Only



Riviera
Glass
Pvt.
Ltd.



C/2116, Oberoi Garden Estates,
Chandivli Studio Road,
Chandivli, Andheri (East),
Mumbai - 400 072, INDIA.
Phone : 2847 5228 / 2847 3297
E-mail : info@rivieraglass.co.in
rivieraglass1995@gmail.com
Website : www.rivieraglass.co.in