



Premium Bottles, Laboratory, with DIN thread, Graduated, with Screw Cap and Pouring Ring (PFA), Thermal Shock Resistant, Autoclavable upto +200 °C, DURAN ISO 4796-1, USP Standard

Premium Bottles contribute to safety in the laboratory while simplifying handling. Important features are

- With enhanced temperature-resistant screw cap manufactured from PFA including a matching pouring ring — can be used from -196 °C to +200 °C. Raw material conforms to USP Class VI — also PTFE coated silicone seal.
- Dye-free cap and pouring ring to prevent marking or contamination of ovens, drying chambers and sterilization cabinets.
- Entire system comprising bottle, cap and pouring ring conforms to USP/FDA.
- Ideal for sterilization processes such as hot-air or dry sterilization and depyrogenation particularly in the pharmaceutical industry.
- Retrace code. • Batch Certificates.



Cat No.	Capacity ml	DIN-thread GL	O. D. mm	Height mm	Qty. Per Pack Pcs
11270 75	100	45	56	100	10
11270 76	250	45	70	138	10
11270 77	500	45	86	176	10
11270 78	1000	45	101	225	10

Bottles, Laboratory, Square, with DIN thread, Graduated, CLEAR, with PP Screw Cap and PP Pouring Ring (Blue), Autoclavable upto 140°C, Suitable for Milk Dilution Application, DURAN USP Standard

Cat No.	Capacity ml	DIN-thread GL	Width mm	Height mm	Qty. Per Pack Pcs
21820 245	100	32	50	109	10
21820 365	250	45	64	143	10
21820 445	500	45	78	181	10
21820 545	1000	45	94	222	10



Bottles, Laboratory, Square, with DIN thread, Graduated, AMBER (Amber processing done in India), with PP Screw Cap and PP Pouring Ring (Blue), Autoclavable upto 140°C, Made from DURAN USP Standard
(Retrace Code Based Batch Certificate shall be provided as per respective sizes of Plain Bottles)

Cat No.	Capacity ml	DIN-thread GL	Width mm	Height mm	Qty. Per Pack Pcs
21826 245	100	32	50	109	10
21826 365	250	45	64	143	10
21826 445	500	45	78	181	10
21826 545	1000	45	94	222	10



Bottles, Laboratory, Square, with DIN thread, Graduated, CLEAR, with Tamper Evident PP Screw Cap, Autoclavable upto 140°C, DURAN USP Standard

Cat No.	Capacity ml	DIN-thread GL	Width mm	Height mm	Qty. Per Pack Pcs
21821 365	250	45	64	143	10
21821 445	500	45	78	181	10
21821 545	1000	45	94	222	10



Bottles, Laboratory, Square, with DIN thread, Graduated, without Cap, Autoclavable upto 140°C, DURAN USP Standard

Cat No.	Capacity ml	DIN-thread GL	Width mm	Height mm	Qty. Per Pack Pcs
21820 24	100	32	50	109	10
21820 36	250	45	64	143	10
21820 44	500	45	78	181	10
21820 54	1000	45	94	222	10

