



DURAN® PURE

DEVELOPED FOR PHARMA

DURAN® PURE is a range of high-quality bottles and closures developed for the needs of the pharmaceutical industry. **DURAN® PURE** meets the requirements of the international pharmacopeias and is available with the certificates required by the regulatory authorities

FEATURES

- **EP/USP Type I borosilicate glass**
- **Size range 25 ml - 20 l**
- **Protective cover** prevents contamination
- **Single batch delivery CoA** available
- **100% traceability**
- **100% visual inspection**
- **Change control**
- **Drug Master File**

APPLICATION AREAS

- Manufacturing of pharmaceutical products
- Storage of bulk drug intermediates
- Transportation of APIs



DURAN PURE bottle with DIN thread GL 45 and dust cover
Cat. No. 81801 5403

Detailed view of dust cover





DURAN® PURE Premium Screw Caps

DURAN® PURE Premium Screw Caps made from PFA (fluoropolymer) are chemically and thermally very resistant and compatible with autoclaving and other common sterilization methods. The products include silicone seals which are coated with highly inert PTFE.



DURAN® PURE Premium Screw Cap GL 45
Cat. No. 89108 8679

FEATURES

- Caps made from PFA with PTFE coated silicone seal
- Produced and packaged in ISO Class 7 clean rooms
- Single batch delivery CoA available
- 100% traceability
- Change control
- Drug Master File



DURAN® PURE Premium Lip Seal Cap GL 45
Cat. No. 89886 2802

Bottles, with DIN thread, Graduated, CLEAR, with Dust Cover without Cap, USP / FDA conformity (DMF No. 28042) DURAN PURE ISO 3585

Cat No.	Capacity ml	DIN-thread GL	O. D. mm	Height mm	Qty Per Pallet Pcs.
81801 1401	25	25	36	74	940
81801 2406	100	45	56	105	540
81801 3608	250	45	70	143	280
81801 4407	500	45	86	181	360
81801 5403	1000	45	101	230	210
81801 6305	2000	45	136	267	120
81801 7301	5000	45	183	340	48
81801 8606	10000	45	230	415	30
81801 9105	20000	45	295	510	12

Bottles, with DIN thread, Graduated, AMBER, with Dust Cover without Cap, USP / FDA conformity (DMF No. 28042) DURAN PURE DIN ISO 3585

Cat No.	Capacity ml	DIN-thread GL	O. D. mm	Height mm	Qty Per Pallet Pcs.
81806 1401	25	25	36	74	940
81806 2406	100	45	56	105	540
81806 3608	250	45	70	143	280
81806 4407	500	45	86	181	360
81806 5403	1000	45	101	230	210
81806 6305	2000	45	136	267	120
81806 7301	5000	45	183	340	48
81806 8606	10000	45	230	415	30
81806 9105	20000	45	295	510	12

Premium Caps Made of PFA, Temperature Resistant from -196 °C to +200 °C, without colourants, USP/FDA conformity, (DMF No. 28044) DURAN PURE

Cat No.	Description	DMF No.	DIN-thread GL	Qty. Per Pack Pcs
89108 1306	Pure Premium Screw Cap	28044	25	5
89108 8679	Pure Premium Screw Cap	28044	45	5
89886 2802	Pure Premium Lip Seal Cap	031063	45	5

DURAN PURE bottles are also available in Wide Mouth GL 80 both in Plain & Amber. Please ask for details.