



Bottles, Laboratory, with DIN thread, Graduated, Plastic Coated, ISO 4796-1, USP Standard with PP Screw Cap and PP Pouring Ring (Blue), Autoclavable upto 140°C, DURAN PROTECT

Cat No.	Capacity ml	DIN-thread GL	O. D. mm	Height mm	Qty. Per Pack Pcs
10926 765	25	25	36	70	10
10926 775	50	32	46	87	10
21805 245	100	45	56	105	10
21805 295	150	45	62	115	10
21805 365	250	45	70	143	10
21805 445	500	45	86	181	10
21805 515	750	45	95	208	10
21805 545	1000	45	101	230	10
21805 635	2000	45	136	267	10
21805 695	3500	45	160	300	1
21805 735	5000	45	182	330	1



Bottles, Laboratory, with DIN thread, Wide Mouth, Graduated, Plastic Coated, USP Standard with PP Screw Cap and PP Pouring Ring (Blue), Autoclavable upto 140°C, DURAN PROTECT

Cat No.	Capacity ml	DIN-thread GL	O. D. mm	Height mm	Qty. Per Pack Pcs
21865 365	250	80	95	105	10
11601 525	500	80	101	148	10
11601 635	1000	80	101	218	10
11601 645	2000	80	136	248	10
21865 695	3500	80	160	271	1
11601 655	5000	80	182	311	1



Bottles, Laboratory, with DIN thread, Graduated, Pressure Resistant, CLEAR, ISO 1595, USP Standard without Screw Cap and Pouring Ring, Autoclavable upto 140°C, DURAN PRESSURE PLUS

Pressure Plus Bottles are the ideal choice where applications involve safe working under pressure or vacuum, sampling under pressure, storage of gas generating media etc. Vacuum and / or pressure resistant from -1 to 1.5 bar due to modified geometry (based on ISO 4796-1). Pressure resistance conforms to EN 1595 (safety tested by TUV). When under full pressure thermal shock resistance is 30K and maximum usage temperature is +140°C. Graduations are printed in blue for visual differentiation with standard laboratory bottles. Retrace code provided on individual bottles for online generation of batch certificates.

Cat No.	Capacity ml	DIN-thread GL	O. D. mm	Height mm	Qty. Per Pack Pcs
21810 24	100	45	56	100	10
10922 34	250	45	70	138	10
10922 35	500	45	86	176	10
21810 54	1000	45	101	225	10



Bottles, Laboratory, with DIN thread, Graduated, Pressure Resistant, AMBER, ISO 1595, USP Standard (Amber processing done in Germany), without Screw Cap and Pouring Ring, Autoclavable upto 140°C, DURAN PRESSURE PLUS

Cat No.	Capacity ml	DIN-thread GL	O. D. mm	Height mm	Qty. Per Pack Pcs
21816 24	100	45	56	100	10
10943 67	250	45	70	138	10
10943 68	500	45	86	176	10
21816 54	1000	45	101	225	10