



Flasks, Volumetric, with Interchangeable Glass Stopper, Accuracy as per **Class 'A', USP,** **White** coloured Marking, with Individual Certificate of Performance, Made from Heat Resistant, Low Expansion Borosilicate Glass, RIVIERA DIN EN ISO 1042



Cat No.	Capacity ml	Tolerance ± ml	Joint Size	Qty. Per Pack Pcs
71022 55	5	0.02	10/19	20
71022 60	10	0.02	10/19	20
71022 65	20	0.03	10/19	20
71022 70	25	0.03	10/19	20
71022 75	50	0.05	10/19	20
71022 80	100	0.08	14/23	20
71022 85	200	0.10	14/23	12
71022 86	250	0.12	14/23	12
71022 90	500	0.15	19/26	12
71022 93	1000	0.30	19/26	12

Flasks, Volumetric, **Class 'A', USP,** with Ergonomic Polyethylene Stopper, **Blue** Graduation, USP Conformity, with USP **Individual Certificate,** DURAN DIN EN ISO 1042



Cat No.	Capacity ml	Tolerance ± ml	Joint Size	Qty. Per Pack Pcs
24671 0958	5	0.02	10/19	2
24671 1054	10	0.02	10/19	2
24671 1457	25	0.03	10/19	2
24671 1757	50	0.05	12/21	2
24671 2556	100	0.08	14/23	2
24671 3252	200	0.10	14/23	2
24671 3655	250	0.12	14/23	2
24671 4454	500	0.20	19/26	2
24671 5459	1000	0.30	24/29	2
24671 6352	2000	0.50	29/32	2

Flasks, Volumetric, **Class 'A', USP,** with Polyethylene Stopper **Blue** Graduation, Conformity Certified, with Batch Certificate - DURAN/EM-TECHCOLOR Brand DIN EN ISO 1042



Cat No.	Capacity ml	Tolerance ± ml	Joint Size	Qty. Per Pack Pcs
28201 5327	5	0.02	10/19	10
28201 6127	10	0.02	10/19	10
28201 6527	20	0.03	10/19	10
28201 7027	25	0.03	10/19	10
28201 7527	50	0.05	12/21	10
28201 8027	100	0.08	12/21	10
28201 8527	200	0.10	14/23	10
28201 8627	250	0.12	14/23	10
28201 8927	500	0.15	19/26	10
28201 9327	1000	0.30	24/29	1
28201 9427	2000	0.50	29/32	1

Available in AMBER - Price on Request

NEW

Riviera Individually calibrated
Volumetric Flasks USP
Volumetric Flasks Class A
Measuring Cylinders Class A
Volumetric Pipettes Class A
Burettes & Automatic Burettes Class A




Now ACCESS your
RIVIERA Volumetric Glassware certification
on your mobile phone...



Riviera Glass - Certificates App



Flasks, Volumetric, with Interchangeable Glass Stopper, Accuracy as per Class 'A', **DIN EN ISO 1042**
Amber Stained Marking, with Individual Certificate of Performance,
 Made from Heat Resistant, Low Expansion Borosilicate Glass, RIVIERA



Cat No.	Capacity ml	Tolerance ±ml	Joint Size	Qty. Per Pack Pcs
71029 55	5	0.025	10/19	20
71029 60	10	0.025	10/19	20
71029 65	20	0.040	10/19	20
71029 70	25	0.040	10/19	20
71029 75	50	0.060	10/19	20
71029 80	100	0.100	14/23	20
71029 85	200	0.150	14/23	12
71029 86	250	0.150	14/23	12
71029 90	500	0.250	19/26	12
71029 93	1000	0.400	19/26	12

Flasks, Volumetric, with Interchangeable Glass Stopper, Accuracy as per Class 'A', **DIN EN ISO 1042**
White coloured Marking, with Individual Certificate of Performance,
 Made from Heat Resistant, Low Expansion Borosilicate Glass, RIVIERA



Cat No.	Capacity ml	Tolerance ±ml	Joint Size	Qty. Per Pack Pcs
71020 55	5	0.025	10/19	20
71020 60	10	0.025	10/19	20
71020 65	20	0.040	10/19	20
71020 70	25	0.040	10/19	20
71020 75	50	0.060	10/19	20
71020 80	100	0.100	14/23	20
71020 85	200	0.150	14/23	12
71020 86	250	0.150	14/23	12
71020 90	500	0.250	19/26	12
71020 93	1000	0.400	19/26	12

Flasks, Volumetric, Class 'A', with Ergonomic Stopper, **Blue** Graduation, **DIN EN ISO 1042**
 with **Individual Certificate** and Certificate of Conformity, DURAN




Cat No.	Capacity ml	Tolerance ±ml	Joint Size	Qty. Per Pack Pcs
24679 0957	5	0.040	10/19	2
24679 1053	10	0.040	10/19	2
24679 1259	20	0.040	10/19	2
24679 1456	25	0.040	10/19	2
24679 1756	50	0.060	12/21	2
24679 2452	100	0.100	12/21	2
24679 3251	200	0.150	14/23	2
24679 3654	250	0.150	14/23	2
24679 4453	500	0.250	19/26	2
24679 5458	1000	0.400	24/29	2
24679 6351	2000	0.600	29/32	2




Flasks, Volumetric, with I/C Glass Stopper, **Amber** Stained Marking, **DIN EN ISO 1042**
 Accuracy as per Class 'B', Made from Heat Resistant,
 Low Expansion Borosilicate Glass, RIVIERA

Cat	Capacity	Tolerance	Joint Size	Qty. Per Pack
No.	ml	± ml		Pcs
71019 55	5	0.05	10/19	20
71019 60	10	0.05	10/19	20
71019 65	20	0.08	10/19	20
71019 70	25	0.08	10/19	20
71019 75	50	0.12	10/19	20
71019 80	100	0.20	14/23	20
71019 85	200	0.30	14/23	12
71019 86	250	0.30	14/23	12
71019 90	500	0.50	19/26	12
71019 93	1000	0.80	19/26	12



Flasks, Volumetric, with I/C Glass Stopper, **White** coloured Marking, **DIN EN ISO 1042**
 Accuracy as per Class 'B', Made from Heat Resistant,
 Low Expansion Borosilicate Glass, RIVIERA

Cat	Capacity	Tolerance	Joint Size	Qty. Per Pack
No.	ml	± ml		Pcs
71000 55	5	0.05	10/19	20
71000 60	10	0.05	10/19	20
71000 65	20	0.08	10/19	20
71000 70	25	0.08	10/19	20
71000 75	50	0.12	10/19	20
71000 80	100	0.20	14/23	20
71000 85	200	0.30	14/23	12
71000 86	250	0.30	14/23	12
71000 90	500	0.50	19/26	12
71000 93	1000	0.80	19/26	12



Flasks, Volumetric, with I/C Glass Stopper, **AMBER**, Accuracy as per Class 'B', **DIN EN ISO 1042**
White coloured Marking, Made from Heat Resistant,
 Low Expansion Borosilicate Glass, RIVIERA

Cat	Capacity	Tolerance	Joint Size	Qty. Per Pack
No.	ml	± ml		Pcs
71010 55	5	0.05	10/19	20
71010 60	10	0.05	10/19	20
71010 65	20	0.08	10/19	20
71010 70	25	0.08	10/19	20
71010 75	50	0.12	10/19	20
71010 80	100	0.20	14/23	20
71010 85	200	0.30	14/23	12
71010 86	250	0.30	14/23	12
71010 90	500	0.50	19/26	12
71010 93	1000	0.80	19/26	12

