

# DURAN® GL 45 BOTTLE CARRYING SYSTEM



The new modular Bottle Carrying System fits standard laboratory bottles that have a GL 45 neck finish and is compatible with four bottle sizes 2, 5, 10 and 20 litres. It is ergonomically designed, and offers health & safety benefits to the user and allows easy and safe transport of bottles around the laboratory or production facility. It makes pouring liquids from the larger bottles more comfortable and safer. All materials are fully autoclavable.

### Carrying Handle

- Carrying handle is suitable for bare or gloved hands
- Handle strap is flexible enough to allow pouring

### Tightening Strap

- With indicator arrow
- End modified to prevent unthreading

### Pouring Handle

- Ergonomically shaped
- Anti-slip grip



### GL 45 Bottle Collar

- The collar is securely attached by pushing down over the GL 45 thread with 4 teeth that lock under transfer ring.
- The collar is only removable by cutting – so attachment is 100% secure in use.
- A cut line indicator for removal of collar when required.

### Horizontal Straps

- Durable stitched construction
- Made from "anti-slip" material

4 sizes available  
2 L / 5 L / 10 L / 20 L



### Fast and secure tightening

of the polypropylene straps with two side strap adjusters for the larger 5, 10 and 20 litre bottles



### Safe and controlled pouring using handle

- Safer pouring of liquids using a second lower "Pouring Handle" located on the side strap
- Available on the 5 litre carriers only



### Easy and safe carrying using handle

- Handles attaches at the neck and fully supports the base



**Lot Specific Retrace Code**  
for traceability and certification

## Bottle Carrying System suitable for GL 45, DURAN

Cat No.	Suitable for Capacity Litres	Qty Per Pack Pcs
29280 6304	2	1
29280 7309	5	1
29280 8605	10	1
29280 9104	20	1