



Cylinders, Measuring, with I/C Joint & Glass Stopper, Class 'B', **DIN EN ISO 4788**
Made from DURAN/EM-TECHCOLOR Brand

Cat	Capacity	Joint Size	Qty. Per Pack
No.	ml		Pcs
223000 60	10	B 10	2
223010 70	25	B 14	2
223010 75	50	B 14	2
223020 80	100	B 19	2
223030 86	250	B 24	2
223050 90	500	B 34	2
223050 93	1000	B 34	1
223050 94	2000	B 34	1



Cylinders, Measuring, with I/C Joint & Glass Stopper, Class 'A', Conformity Certified, with **DIN EN ISO 4788**
Batch Certificate, Calibrated at 3 Graduation levels for accuracy and tolerance, Made from DURAN/EM-TECHCOLOR Brand
 Batch certificate shall be for Measuring Cylinder Class 'A' only (Cat. No. 22401..)

Cat	Capacity	Joint Size	Qty. Per Pack
No.	ml		Pcs
224000 60	10	B 10	2
224010 70	25	B 14	2
224010 75	50	B 14	2
224020 80	100	B 19	2
224030 86	250	B 24	2
224050 90	500	B 34	2
224050 93	1000	B 34	1
224050 94	2000	B 34	1

Desiccators, Plain, Clear Glass, RIVIERA

Cat	Max. I. D.	Overall Height	Qty. Per Pack
No.	mm	mm	Pcs
74300 46	120	180	1
74300 57	160	230	1
74300 61	210	300	1
74300 66	240	340	1
74300 69	300	420	1



Desiccators, with Plane Flanges, without Junction, with Knobbed Lid, DURAN

Cat	DN	Height	Qty. Per Pack
No.		mm	Pcs
24781 46	100	187	1
24781 57	150	252	1
24781 61	200	309	1
24781 66	250	357	1
24781 69	300	433	1




Desiccators, Vacuum, Clear Glass, RIVIERA

Cat	Max. I. D.	Overall Height	Qty. Per Pack
No.	mm	mm	Pcs
74310 46	120	227	1
74310 57	160	285	1
74310 61	210	347	1
74310 66	240	381	1
74310 69	300	467	1





Desiccators, Vacuum, with Flat Flanges and NOVUS NS Junction, Outlet in Lid, with Stopcock, DURAN				DIN 12491
Cat	DN	Height	Qty. Per Pack	
No.		mm	Pcs	
24782 46	100	174	1	
24782 57	150	239	1	
24782 61	200	296	1	
24782 66	250	344	1	
24782 69	300	420	1	



Porcelain Plates for All Desiccators, RIVIERA			
Cat	DN	Qty. Per Pack	
No.		Pcs	
74320 46	100	1	
74320 57	150	1	
74320 61	200	1	
74320 66	240	1	
74320 69	300	1	



Base for Plane and Vacuum Desiccator, DURAN				DIN 12491
Cat	DN	Height	Qty. Per Pack	
No.		mm	Pcs	
24770 46	100	58	1	
24770 57	150	81	1	
24770 61	200	115	1	
24770 66	250	120	1	
24770 69	300	150	1	



Lid for Plane Desiccator, DURAN				
Cat	DN	Height	Qty. Per Pack	
No.		mm	Pcs	
24410 46	100	75	1	
24410 57	150	98	1	
24410 61	200	107	1	
24410 66	250	122	1	
24410 69	300	150	1	



Lid for Vacuum Desiccator, DURAN				DIN 12491
Cat	DN	Joint Size	Height	Qty. Per Pack
No.			mm	Pcs
24420 46	100	B 24	62	1
24420 57	150	B 24	85	1
24420 61	200	B 24	94	1
24420 66	250	B 24	109	1
24420 69	300	B 24	137	1

Stopcock with PTFE Spindle for Vacuum Desiccator			
Cat	Joint Size	Length	Qty. Per Pack
No.		mm	Pcs
24799 04	B 24	124	1



Dishes, Crystallizing, Made from Heat Resistant, Low Expansion Borosilicate Glass, RIVIERA			
Cat	Capacity	O. D. x Height	Qty. Per Pack
No.	ml	mm	Pcs
71600 44	300	95 x 55	6
71600 54	900	140 x 75	4
71600 59	2000	190 x 90	1



Dishes, Crystallizing, with Spout, DURAN			DIN 12338
Cat	Capacity	O. D. x Height	Qty. Per Pack
No.	ml	mm	Pcs
21311 44	300	95 x 55	10
21311 54	900	140 x 75	10
21311 59	2000	190 x 90	10



Dishes, Evaporating, Made from Heat Resistant, Low Expansion Borosilicate Glass, RIVIERA			
Cat	Capacity	Max. O. D. x Height	Qty. Per Pack
No.	ml	mm	Pcs
71700 41	90	80 x 45	10
71700 44	170	105 x 55	2
71700 54	600	150 x 80	2



Dishes, Evaporating, with Spout, DURAN			DIN 12336
Cat	Capacity	Max. O. D. x Height	Qty. Per Pack
No.	ml	mm	Pcs
21301 41	90	80 x 45	10
21301 44	170	95 x 55	10
21301 54	600	140 x 80	10



Dishes, Petri, Made from Heat Resistant, Low Expansion Borosilicate Glass, RIVIERA		
Cat	O. D. x Height	Qty. Per Pack
No.	mm	Pcs
71500 40	50 x 17	10
71500 42	80 x 17	10
71500 48	100 x 17	10
71500 52	150 x 20	10
71500 59	200 x 20	10

