



**Bottles, Weighing, with I/C Stopper,
Made from Heat Resistant, Low Expansion Borosilicate Glass, RIVIERA**

Cat	Capacity	App. O. D. x Ht. of Bottle only	Qty. Per Pack
No.	ml	mm	Pcs
70900 05	5	20 x 40	10
70900 15	15	25 x 65	10
70900 25	25	30 x 65	10
70900 60	60	40 x 90	4




**Bottles, Weighing, Squat Form, with I/C Stopper,
Made from Heat Resistant, Low Expansion Borosilicate Glass, RIVIERA**

Cat	Capacity	App. O. D. x Ht. of Bottle only	Qty. Per Pack
No.	ml	mm	Pcs
70900 20	20	50 x 35	10
70900 40	40	60 x 40	10




**Bottles, Wash, Complete with I/C Pressure and Delivery Tube,
Made from Heat Resistant, Low Expansion Borosilicate Glass, RIVIERA**

Cat	Capacity	Joint size	Qty. Per Pack
No.	ml		Pcs
74210 36	250	24/29	10
74210 44	500	24/29	10
74210 54	1000	29/32	10



**I/C Pressure and Delivery Tube only for Wash Bottle,
Made from Heat Resistant, Low Expansion Borosilicate Glass, RIVIERA**

Cat	Capacity	Joint size	Qty. Per Pack
No.	ml		Pcs
74211 36	250	24/29	10
74211 44	500	24/29	10
74211 54	1000	29/32	10



**Bottles, Gas Washing, Complete with I/C Pressure and Delivery Tube,
Made from Heat Resistant, Low Expansion Borosilicate Glass, RIVIERA**

Cat	Capacity	Joint size	Qty. Per Pack
No.	ml		Pcs
74212 25	125	29/32	10
74212 36	250	29/32	10
74212 44	500	29/32	10



**I/C Pressure and Delivery Tube only for Gas Washing Bottle,
Made from Heat Resistant, Low Expansion Borosilicate Glass, RIVIERA**

Cat	Capacity	Joint size	Qty. Per Pack
No.	ml		Pcs
74213 25	125	29/32	10
74213 36	250	29/32	10
74213 44	500	29/32	10

