



**Bottles, Aspirator, Graduated, with GL 45 Screwthread and Autoclavable Cap Complete with NS Stopcock and Joint, Made from DURAN Laboratory Bottles**



Cat No.	Capacity ml	Bottle O. D. mm	Total Height mm	Qty. Per Pack Pcs
24703 54	1000	101	230	1
24703 63	2000	136	267	1
24703 73	5000	183	340	1
24703 86	10000	230	415	1
24703 91	20000	295	510	1

**Bottles, B.O.D., with I/C Stopper, RIVIERA**

Suitable for determination of B.O.D. in industrial effluents, wastes and water as per method specified by the American Public Health Association



Cat No.	Capacity ml	App. O.D.x Height of Bottle only mm	Qty. Per Pack Pcs
74100 365	300	73 x 131	10

**Bottles Reagent, Narrow Mouth, Plain, Made from Heat Resistant, Low Expansion Borosilicate Glass, RIVIERA**



Cat No.	Capacity ml	Qty. Per Pack Pcs
70200 18	60	10
70200 25	125	10
70200 36	250	10
70200 44	500	10
70200 54	1000	8
70200 63	2000	8
70200 65	2500	1
70200 73	5000	1
70200 86	10000	1
70200 91	20000	1

**Bottles Reagent, Narrow Mouth, AMBER, Made from Heat Resistant, Low Expansion Borosilicate Glass, RIVIERA**



Cat No.	Capacity ml	Qty. Per Pack Pcs
70210 18	60	10
70210 25	125	10
70210 36	250	10
70210 44	500	10
70210 54	1000	8
70210 65	2500	1