

YOUTILITY BOTTLES LABELS

YOUTILITY LABELS were developed to address the need for a reliable, short term method of sample and bottle tracking and identification. The labels will stay on when you need them to do so, and are easy to remove when no longer required. This allows the bottles to be quickly and easily cleaned, ready for re-use.

The **YOUTILITY** label, is large rectangular label that fits all four sizes of bottle. The labeling of bottles helps to prevent sample mix-ups.

YOUTILITY label formats are manufactured from a tear resistant, matt white, polyester material coated with a high performance adhesive. The print receptive top surface will accept a variety of inks, including copier, printer, ball point and

technical markers. The labels are chemically resistant to typical laboratory chemicals such as disinfectants and many solvents. They have a wide -40 to + 150°C temperature performance range making them suitable for use in fridges, freezers, autoclaves, incubators and water baths without the risk of the labels drying out or falling off.


The **YOUTILITY** labels conveniently combine the flexibility of a blank written label with an autoclave/steam sterilization Class 1 process indicator (ANSI/AAMI/ISO 11140-1:2005). The labels are non-hazardous, and free of toxic heavy metals.

After use, the **YOUTILITY** labels may be easily and completely removed by peeling. The labels are supplied in a re-sealable, splash resistant box containing 100 labels.



Bottles, YOUTILITY, Clear Glass, with Screw Cap, Pouring Ring and Bottle Tag, Autoclavable upto 140°C, DURAN						USP Standard
Cat No.	Capacity ml	DIN-thread GL	O. D. mm	Height mm	Qty. Per Pack Pcs	
21881 285	125	45	55	124	4	
21881 365	250	45	66	158	4	
21881 445	500	45	78	193	4	
21881 545	1000	45	93	253	4	

Utility bottles, AMBER colour processing done either in Germany or in India - Prices available On Request



Screw Cap & Pouring Ring, PP, Cyan, suitable for YOUTILITY Bottles, DURAN			
Cat No.	Description	DIN Thread GL	Qty. Per Pack Pcs
29229 28	Cap	Food Grade	10
29241 28	Pouring Ring	—	16



Tags in Silicone for YOUTILITY Bottles, DURAN			
Cat No.	Colour	Qty. Per Pack Pcs	
29243 29	Multi colour	16	
29243 2818	Red	20	
29243 2826	Orange	20	
29243 2834	Yellow	20	
29243 2842	Green	20	
29243 2859	Blue	20	
29243 2867	Purple	20	
29243 2875	Black	20	
29243 2883	White	20	



Self Adhesive Label (Bottles only), for YOUTILITY Bottles, DURAN	
Cat No.	Qty. Sets
29401 0203	100 Sets