



Tubes, Centrifuge, Conical, Graduated with White Graduation, HIRSCHMANN

Cat	Capacity	O. D.	Total Length ± 5 mm	Qty. Per Pack Pcs
No.	ml	mm		
60201 01	15	0.1	112	100

**Test Tubes, with Rim, Borosilicate Glass,
Made from Heat Resistant, Low Expansion Borosilicate Glass, (suitable for heating applications), RIVIERA**

Cat	Size	Qty. Per Pack Pcs
No.	mm	
71100 001	10 x 75	100
71100 002	12 x 75	100
71100 003	12 x 100	100
71100 004	15 x 125	100
71100 005	15 x 150	100
71100 006	18 x 150	100
71100 007	25 x 100	50
71100 008	25 x 150	50
71100 009	25 x 200	50
71100 010	32 x 200	25
71100 011	38 x 200	25

**Test Tubes, without Rim, Borosilicate Glass,
Made from Heat Resistant, Low Expansion Borosilicate Glass, (suitable for heating applications), RIVIERA**

Cat	Size	Qty. Per Pack Pcs
No.	mm	
71110 001	10 x 75	100
71110 002	12 x 75	100
71110 003	12 x 100	100
71110 004	15 x 125	100
71110 005	15 x 150	100
71110 006	18 x 150	100
71110 007	25 x 100	50
71110 008	25 x 150	50
71110 009	25 x 200	50
71110 010	32 x 200	25
71110 011	38 x 200	25

Test Tubes, with Rim, Borosilicate Glass, (NOT suitable for heating applications), RIVIERA

Cat	Size	Qty. Per Pack Pcs
No.	mm	
23277 001	10 x 75	100
23277 002	12 x 75	100
23277 003	12 x 100	100

Test Tubes, without Rim, Borosilicate Glass, (NOT suitable for heating applications), RIVIERA

Cat	Size	Qty. Per Pack Pcs
No.	mm	
23278 001	10 x 75	100
23278 002	12 x 75	100
23278 003	12 x 100	100

