

Bottle Top Liquid Dispenser, with Ceramic Piston, Fixed Volume, Autoclavable at 121° C, 2 Bar, Conformity Certified, with Individual Number & Work Certificate, CERAMUS - CLASSIC FIX, make HIRSCHMANN.

Equipment shall be supplied complete with Thread A 32 mm, 3 thread adapters (A 28, S 40 & A 45), 1 screw coupled suction tube and 1 discharge tube unit complete

Cat No.	Volume ml	Qty. Per Pack Pcs
93111 00	1	1
93211 00	2	1
93311 00	5	1
93411 00	10	1
93511 00	25	1
93611 00	50	1



Bottle Top Liquid Dispenser, with Ceramic Piston, Variable Volume, Autoclavable at 121° C, 2 Bar, Conformity Certified, with Individual Number & Work Certificate, CERAMUS - CLASSIC, make HIRSCHMANN.

Equipment shall be supplied complete with Thread A 32 mm, 3 thread adapters (A 28, S 40 & A 45), 1 screw coupled suction tube and 1 discharge tube unit complete.

Cat No.	Volume ml	Qty. Per Pack Pcs
93100 00	0.2 - 1	1
93200 00	0.4 - 2	1
93300 00	1.0 - 5	1
93400 00	2.0 - 10	1
93500 00	5.0 - 30	1
93600 00	10.0 - 60	1



Bottle Top Liquid Dispenser, Variable Volume, with Ceramic Piston & Reagent recirculation system, Autoclavable at 121° C, 2 Bar, Conformity Certified, with Individual Number & Work Certificate, CERAMUS®, make HIRSCHMANN

Equipment shall be supplied suitable for A 45 mm Thread, 3 thread adapters (A 32, S 40 & A 38), 1 screw coupled suction tube, recirculation tube, spanner and 1 discharge tube unit complete.

Cat No.	Volume ml	Volume Division ml	Accuracy R%	Reproducibility CV%	Qty. Per Pack Pcs
93120 00	0.2 - 1	0.05	0.6	0.2	1
93220 00	0.4 - 2	0.05	0.5	0.1	1
93320 00	1 - 5	0.10	0.5	0.1	1
93420 00	2 - 10	0.25	0.5	0.1	1
93520 00	5 - 30	0.50	0.5	0.1	1
93620 00	10 - 60	1.00	0.5	0.1	1

ceramus®
with unrivalled
ceramic piston



ceramus® classic
dispensing as perfect
as it could be

