

# Bottle Top Dispenser - RIVIERA



Without Recirculation



With Recirculation

## Features :

- Ergonomic Design enables low pressure dispensing for even larger volumes
- Dual graduation marks enable precise volume setting & visibility from both sides
- Unique piston mechanism allows cleaning & servicing of piston cylinder without disturbing calibration
- Volume lock ensures any accidental discharge
- Spring-less valve design offers high-chemical resistance and long-life with minimal maintenance
- Telescopic filling tube suitable for different sizes of Bottles
- Discharge tube with 360° rotation makes it convenient to dispense in multiple vessels
- The device comes with 4 additional adapters for most common bottle sizes (28mm, 38mm, 40mm, 45mm)
- The device can be easily calibrated with a single tool
- Calibrated as per ISO 8655 standard. Provided with individual calibration certificate

### Bottle Top Dispenser, without Recirculation, RIVIERA

Cat No.	Recirculation	Volume Range ml	Volume Division ml	Accuracy %	Repeat-ability %	Qty Pack Per Pcs.
78600 0025	No	0.25 - 2.5	0.05	0.6	0.2	1
78600 0050	No	0.5 - 5	0.10	0.5	0.2	1
78600 0100	No	1- 10	0.20	0.5	0.2	1
78600 0250	No	2.5 - 25	0.50	0.5	0.2	1
78600 0300	No	2.5 - 30	0.50	0.5	0.2	1
78600 0500	No	5 - 50	1.00	0.5	0.2	1
78600 0600	No	5 - 60	1.00	0.5	0.2	1

### Bottle Top Dispenser, with Recirculation RIVIERA

Cat No.	Recirculation	Volume Range ml	Volume Division ml	Accuracy %	Repeat-ability %	Qty Pack Per Pcs.
78700 0025	Yes	0.25 - 2.5	0.05	0.6	0.2	1
78700 0050	Yes	0.5 - 5	0.10	0.5	0.2	1
78700 0100	Yes	1- 10	0.20	0.5	0.2	1
78700 0250	Yes	2.5 - 25	0.50	0.5	0.2	1
78700 0300	Yes	2.5 - 30	0.50	0.5	0.2	1
78700 0500	Yes	5 - 50	1.00	0.5	0.2	1
78700 0600	Yes	5 - 60	1.00	0.5	0.2	1