



**MICROPIPETTES, Fixed Volume, Single Channel, Fully Autoclavable, with Individual Work Certificate, RIVIERA +**

Does not require recalibration after Autoclaving. Fully Autoclavable since material of construction of body is high grade engineering plastics. Colour coded for easy identification. Calibration conforming to DIN 12650 & EN-ISO 8655 standards ensuring high accuracy & precision. Calibration report with every pipette.

Cat No.	Volume $\mu$ l	Qty. Per Pack Pcs
77200 01	5	1
77200 02	10	1
77200 03	20	1
77200 04	25	1
77200 05	50	1
77200 06	100	1
77200 07	200	1
77200 08	250	1
77200 09	500	1
77200 10	1000	1
77200 11	2000	1
77200 12	5000	1
77200 13	10000	1



**MICROPIPETTES, Variable Volume, Single Channel, Fully Autoclavable, with Individual Work Certificate, RIVIERA +**

Does not require recalibration after Autoclaving. Fully Autoclavable since material of construction of body is high grade engineering plastics. Colour coded for easy identification. Calibration conforming to DIN 12650 & EN-ISO 8655 standards ensuring high accuracy & precision. Calibration report with every pipette.

Cat No.	Volume $\mu$ l	Volume Division $\mu$ l	Qty. Per Pack Pcs
77300 01	0.5 - 10	0.1	1
77300 02	5 - 50	0.5	1
77300 03	50 - 200	1.0	1
77300 04	200 - 1000	5.0	1
77300 05	10 - 100	0.5	1
77300 06	100 - 1000	5.0	1
77300 07	2 - 20	0.1	1
77300 08	20 - 200	1.0	1
77300 09	1000 - 5000	50.0	1
77300 10	1000 - 10000	100.0	1

**Plastic Micro-Pipette Tips, Universal Type**



Cat No.	Capacity Range $\mu$ l	Compatible with Micro-Pipette	Qty per Pack
77400 001	0.5 - 20	77301 01, 77301 07, 77300 01, 77300 07	1000
77400 002	0.5 - 200	77200 01, 77200 02, 77200 03, 77200 04, 77200 05, 77200 06, 77200 07, 77201 01, 77201 02, 77201 03, 77201 04, 77201 05, 77201 06, 77201 07, 77300 02, 77300 03, 77300 05, 77300 08, 77301 02, 77301 03, 77301 05, 77301 08	1000
77400 003	100 - 1000	77200 08, 77200 09, 77200 10, 77300 04, 77300 06, 77201 08, 77201 09, 77201 10, 77301 04, 77301 06	1000
77400 004	1000 - 5000	77200 11, 77200 12, 77300 09, 77201 11, 77201 12, 77301 09	100
77400 005	1000 - 10000	77201 13, 77301 10, 77200 13, 77300 10	100