



MICROPIPETTES, Fixed Volume, Single Channel, with Individual Work Certificate, RIVIERA

Recalibration required after each Autoclaving. Colour coded for easy identification. Calibration conforming to DIN 12650 & EN-ISO 8655 standards ensuring high accuracy & precision. Calibration report with every pipette.

Cat No.	Volume μ l	Qty. Per Pack Pcs
77201 01	5	1
77201 02	10	1
77201 03	20	1
77201 04	25	1
77201 05	50	1
77201 06	100	1
77201 07	200	1
77201 08	250	1
77201 09	500	1
77201 10	1000	1
77201 11	2000	1
77201 12	5000	1
77201 13	10000	1



MICROPIPETTES, Variable Volume, Single Channel, with Individual Work Certificate, RIVIERA

Recalibration required after each Autoclaving. Colour coded for easy identification. Calibration conforming to DIN 12650 & EN-ISO 8655 standards ensuring high accuracy & precision. Calibration report with every pipette.

Cat No.	Volume μ l	Volume Division μ l	Qty. Per Pack Pcs
77301 01	0.5 - 10	0.1	1
77301 02	5 - 50	0.5	1
77301 03	50 - 200	1.0	1
77301 04	200 - 1000	5.0	1
77301 05	10 - 100	0.5	1
77301 06	100 - 1000	5.0	1
77301 07	2 - 20	0.1	1
77301 08	20 - 200	1.0	1
77301 09	1000 - 5000	50.0	1

Plastic Stand for Micro Pipettors

Cat No.	To Hold No. of Micro Pipettor	Qty. Per Pack Pcs
77500 001	4	1
77500 002	6	1

