

Pipette Controller - RIVIERA



Main Features

- ✓ Easy to use Operating Knobs
- ✓ Powerful Pump
- ✓ Double Safety Valve and Hydrophobic Filters
 - Just 4 seconds to fill a 25ml Pipette
 - No fluid penetration ensured by Double Safety Valve and Hydrophobic Filters
 - Fast and Slow Speed modes
 - Excellent Control of Speed through finger pressure control
 - Low Battery Indicator
 - Power through Cadmium free Batteries - Environment friendly – 2 x 1.2V/1000mAh NiMH batteries. Can be used continuously for 8 hrs. Battery Charging 2-3 hrs.
 - Range of Pipette volume from 0.1 ml to 200ml

Accessories available with the Controller

- ✓ Battery Charger
- ✓ Hydrophobic PTFE membrane filter 0.45 μm 1pc
- ✓ Hydrophobic PTFE membrane filter 0.20 μm 1pc
- ✓ Wall mounting stand with screw

Pipette Controller, RIVIERA

Cat No.	Equipment Type	Operating Principle	Operating Range ml	Accessories Supplied	Qty. Per Pack Pcs
77510 000	Pipette Controller	Electrical	0.1 - 200	As above	1