

## Features :

- No Moving parts -truly maintenance free
- Low aspect ratio (just 13/14 mm thick)
- Software driven gradual acceleration
- Microprocessor controlled with last run memory feature
- Pulse mode-reverses the direction every 30 seconds for proper mixing
- Universal power supply –small wall mount power adapter
- IP 65 compliant design (water spill resistant)
- Corrosion resistant SS top surface to fight harsh chemicals
- Synchronized speed setting (Multi-position)



Specifications	iStir UNO	iStir Quattro
Max RPM	1500	1200
Max Setting	From 15 RPM to every 50 RPM	
Capacity ml (H <sub>2</sub> O)	1 X 800	4 x 400
Time setting	1 to 99 min	
Digital Display	Yes	
Recommended bar length mm	25	
Pulse function	Yes	
Sync Mode	NA	Yes
Lat run memory	Yes	
Dimensions (W x D x H) in mm	124 x 191 x 14.5	457 x 126 x 13.4
Power Consumption	5 W	13 W
Weight	400 gm	1.5 Kg
Protection class	IP 65 compliant design (water spill resistant)	

Cat No.	Product Code	Description
3006 100	MS01-1111	Motor less ultra slim single position magnetic stirrer with 1500 RPM
3006 200	MS01-4110	Motor less ultra slim 4 position magnetic stirrer with 1200 RPM

### Optional Accessories

Cat No.	Product Code	Description
3006 000	SB25	Magnetic bar (L X D) 25x 8 mm