

Features : (Mini Vortex Mixers)

- Tiny Vortex with small footprint
- BLDC motor for maintenance free long life
- Smart & robust design- Press activation & suspension
- Robust construction with dia cast aluminum base and rubber pads for stable & salient operation



Features : (Vortex Mixers)

- Brushless DC motor fitted with ball bearing for efficient noiseless operation
- Adjustable speed range with accuracy of $\pm 5\%$
- Time setting 0 to 999 min (iSwix VT)
- Programmable pulse mode function (iSwix VT)
- 3-way switch for touch, continuous & stand-by mode
- Visual indication marker for max speed

Specifications	iSwix Jr	iSwix Jr VT	iSwix	iSwix VT
Max RPM	4500		2800	4200
Orbital diameter mm	4.5		4.0	
Motor type	BLDC		AC Synchronous	BLDC
Digital Display	No		No	Yes
Max Load Capacity gm	100		500	
Max tube diameter	30 mm test tube		Optional Tubes adapters available	
Press activation	Yes		3 way switch	
Speed setting	No	Yes	Yes	
Time setting	No		No	Yes
Dimensions (W x D x H) in mm	100 x 100 x 70		205 x 136 x 138.5	
Power supply	24 VDC		220 VAC, 50 Hz	24 VDC
Power Consumption	14 W		20 W	
Weight	0.5 kg		3 kg	

Cat No.	Product Code	Description
3003 100	iSwix Jr	Analog tiny vortex mixer with fixed speed
3004 100	iSwix Jr VT	Analog tiny vortex mixer with variable speed
3003 200	iSwix	Analog vortex mixer with variable speed, 230 VAC operated
3004 200	iSwix VT	Digital vortex mixer with variable speed & time

Optional Accessories

Cat No.	Product Code	Description
3005 001	VM-A1	Universal Attachment
3005 002	VM-A3	Micro plate attachment for iSwix/iSwix VT Vortex Mixer
3005 003	VM-A4	Flask attachment for iSwix/iSwix VT Vortex Mixer
3005 004	VM-A2	Micro Centrifuge Tubes attachment for iSwix/iSwix VT Vortex Mixer
3005 005	VM-A6	Round Disc for for iSwix/iSwix VT Vortex Mixer