



**Distillation Unit, Single Stage, with Borosilicate Boiler & Condenser, Built-in Cut-Off with Reset Button, CABINET TYPE, Horizontal Model, RIVIERA**

Cat No.	Type	Power kW	Cooling Water Requirement Ltrs/Hr (Approx.)	Output Ltrs. / Hr	Conductivity $\mu$ s/cm	Qty. Per Pack Pcs
72212 15	Single	1.5	35	1.5	1.5 - 2.5	1
72201 157	Condenser only for above unit					
72201 1555	Boiler only for above unit with PT 100 Sensor provision					
72201 1566	Quartz Heater only for above unit					
72212 40	Single	3.0	50	4.0	1.5 - 2.5	1
722210 043	Condenser only for above unit					
722210 0415	Boiler only for above unit with PT 100 Sensor provision					
722210 0426	Ceramic Heater only for above unit					

**Distillation Unit, Double Stage, with Borosilicate Boiler & Condenser, Built-in Cut-Off with Reset Button, CABINET TYPE, Horizontal Model, RIVIERA**

Cat No.	Type	Power kW	Cooling Water Requirement Ltrs/Hr (Approx.)	Output Ltrs. / Hr	Conductivity $\mu$ s/cm	Qty. Per Pack Pcs
72213 15	Double	3.0	35 x 2 inputs	1.5	0.8 - 1.2	1
72201 157	Condenser only for above unit					
72201 1555	Boiler only for above unit with PT 100 Sensor provision					
72201 1566	Quartz Heater only for above unit					
72213 40	Double	6.0	50 x 2 inputs	4.0	0.8 - 1.2	1
722210 043	Condenser only for above unit					
722210 0415	Boiler only for above unit with PT 100 Sensor provision					
722210 0426	Ceramic Heater only for above unit					

**Distillation Unit, Single Stage, All Quartz, with Built-in Quartz Heater, Vertical Panel Mounted, RIVIERA**

Cat No.	Type	Power kW	Cooling Water Requirement Ltrs/Hr (Approx.)	Output Ltrs. / Hr	Conductivity $\mu$ s/cm	Qty. Per Pack Pcs
72220 10	Single	1.0	30	1.0	1.5 - 2.0	1
72220 105	Boiler only for above unit					
72220 106	Condenser only for above unit					
72220 20	Single	1.5	40	2.0	1.5 - 2.0	1
72220 205	Boiler only for above unit					
72220 206	Condenser only for above unit					
72220 30	Single	2.0	50	3.0	1.5 - 2.0	1
72220 305	Boiler only for above unit					
72220 306	Condenser only for above unit					
72220 40	Single	3.0	60	4.0	1.5 - 2.0	1
72220 405	Boiler only for above unit					
72220 406	Condenser only for above unit					
72220 50	Single	4.0	65	5.0	1.5 - 2.0	1
72220 505	Boiler only for above unit					
72220 506	Condenser only for above unit					

**Note : The above quality of Distillate is obtained by using treated water as boiler feed & collecting the Distillate in closed quartz container.**

**Riviera** Water Distillation Units



**Cat. No. 72240 Series**



**Cat. No. 72260 000**



**Cat. No. 72230 Series**



**Cat. No. 72211 Series**



**Cat. No. 72212 Series**