




Distillation Unit, All Glass, with heater and In- Built Cut-off System, Horizontal Type, Made from Heat Resistant, Low Expansion Borosilicate Glas



Specifications:

- Single stage/Double Stage
- Horizontal model
- Boiler & Condenser made of Heat Resistant, low expansion Borosilicate Glass with Replaceable Quartz/Ceramic heater.
- Output Capacity: 1.5 & 4.0 Ltrs/hr
- Electrical Requirement: 230 VAC, 50 Hz, Single Phase
- Biological Activity: Pyrogen Free
- Conductivity: 1.5 to 2.5 $\mu\text{S/cm}$
- Temperature of Distillate: 35°C to 40°C

Ordering Information:

Cat. No.	Stages	Type of Boiler	Type of Condenser	Cooling water requirement Ltrs/Hr. (Approx.)	Output Ltrs/Hr.	Conductivity $\mu\text{S/cm}$	Power kW
72214 15	Single	Borosilicate	Borosilicate	35	1.5	1.5 – 2.5	1.25
72214 40	Single	Borosilicate	Borosilicate	50	4.0	1.5 – 2.5	3.0
72214 1515	Double	Borosilicate	Borosilicate	35 x 2 inputs	1.5	0.8 – 1.2	3.0
72214 4040	Double	Borosilicate	Borosilicate	50 x 2 inputs	4.0	0.8 – 1.2	6.0

=====

Riviera Glass Pvt. Ltd. | C-2116, Oberoi Garden Estates, Chandivali, Andheri (E), Mumbai – 400072 Tel.No. 022-28475228 | Email- info@rivieraglass.co.in, rivieraglass1995@gmail.com | www.rivieraglass.co.in