

**Riviera** Water Distillation Units



**Cat. No. 72240 Series**



**Cat. No. 72260 000**



**Cat. No. 72230 Series**



**Cat. No. 72211 Series**



**Cat. No. 72212 Series**



**Distillation Unit, Single Stage, with Quartz Boiler, Borosilicate Condenser & Built-in Quartz Heater,  
Vertical Panel Mounted, RIVIERA**

Cat No.	Type	Power kW	Cooling Water Requirement Ltrs/Hr (Approx.)	Output Ltrs. / Hr	Conductivity $\mu$ s/cm	Qty. Per Pack Pcs
72230 010	Single	1.0	30	1.0	1.5 - 2.0	1
72230 110	Boiler, Quartz only for above unit					1
72230 210	Condenser, Borosilicate only for above unit					1
72230 020	Single	1.5	40	2.0	1.5 - 2.0	1
72230 120	Boiler, Quartz only for above unit					1
72230 220	Condenser, Borosilicate only for above unit					1
72230 030	Single	2.0	50	3.0	1.5 - 2.0	1
72230 130	Boiler, Quartz only for above unit					1
72230 230	Condenser, Borosilicate only for above unit					1
72230 040	Single	3.0	60	4.0	1.5 - 2.0	1
72230 140	Boiler, Quartz only for above unit					1
72230 240	Condenser, Borosilicate only for above unit					1
72230 050	Single	4.0	65	5.0	1.5 - 2.0	1
72230 150	Boiler, Quartz only for above unit					1
72230 250	Condenser, Borosilicate only for above unit					1

**Distillation Unit, Double Stage, with Quartz Boiler, Quartz Condenser & Built-in Quartz Heater,  
Vertical Panel Mounted, RIVIERA**

Cat No.	Type	Power kW	Cooling Water Requirement Ltrs/Hr (Approx.)	Output Ltrs. / Hr	Conductivity $\mu$ s/cm	Qty. Per Pack Pcs
72240 010	Double	2.0	40	1.0	1.5 - 2.0	1
72240 011	Lower Boiler, Quartz only for above unit					1
72240 012	Upper Boiler with Condenser, Quartz only for above unit					1
72240 015	Double	3.0	60	1.5	1.5 - 2.0	1
72240 0151	Lower Boiler, Quartz only for above unit					1
72240 0152	Upper Boiler with Condenser, Quartz only for above unit					1
72240 020	Double	4.0	60	2.0	1.5 - 2.0	1
72240 021	Lower Boiler, Quartz only for above unit					1
72240 022	Upper Boiler with Condenser, Quartz only for above unit					1
72240 025	Double	4.5	70	2.5	1.5 - 2.0	1
72240 0251	Lower Boiler, Quartz only for above unit					1
72240 0252	Upper Boiler with Condenser, Quartz only for above unit					1
72240 030	Double	5.0	75	3.0	1.5 - 2.0	1
72240 031	Lower Boiler, Quartz only for above unit					1
72240 032	Upper Boiler with Condenser, Quartz only for above unit					1
72240 040	Double	6.0	85	4.0	1.5 - 2.0	1
72240 041	Lower Boiler, Quartz only for above unit					1
72240 042	Upper Boiler with Condenser, Quartz only for above unit					1
72240 050	Double	8.0	100	5.0	1.5 - 2.0	1
72240 051	Lower Boiler, Quartz only for above unit					1
72240 052	Upper Boiler with Condenser, Quartz only for above unit					1

**Note : The above quality of Distillate is obtained by using treated water as boiler feed & collecting the Distillate in closed quartz container.**