

Precision Balance CG Series



- Auto calibration with external weight.
- Easy to read LCD Display.
- Multifunction weighing units.
- Auto zero tracking, low battery for continues use.
- Overload protection design.
- S.S. Weighing Pan level indicator Quick weighing accurate result.

Model	CG 203L	CG 303L
Capacity	200 gm	300 gm
Readability	0.001 gm	0.001 gm
Repeatability (+/-)	0.001 gm	0.001 gm
Linearity (+/-)	0.002 gm	0.002 gm
Pan Size (mm/inc)	90 Ø	90 Ø
Cat No.	1067 0203	1067 0303



Standard Features

- Compact and space saving.
- Large Stainless steel plate.
- Accurate to within \pm division at full load.
- Overload Protection.
- Dart Shield for more accurate reading.

- Standard RS 232 Interface for data output.
- LED 7 Digit Display for easy viewing.
- Standard dual side or External Display.
- Auto-Zero Tracking, Automatic External Calibration
- 8 Different weighing units such as g, ct, Tola etc.
- Piece Counting & Percentage Weighing Function
- With leveling bubble and adjustable feet for the stability of scale.
- Inbuilt rechargeable battery for continues use.
- AC or DC Power supply : Output AC/DC Adapter.

Model	CG 302	CG 602	CG 1202
Capacity	300 gm	600 gm	1200 gm
Readability	0.01 gm	0.01 gm	0.01 gm
Repeatability (sd)	\pm 0.03 gm	\pm 0.03 gm	\pm 0.03 gm
Linearity (+/-)	0.01 gm	0.01 gm	0.01 gm
Pan Size (mm)	118 Ø		145 x 125
Cat No.	1066 0302	1066 0602	1066 1202