



KERN ECB • EOB

Platform scales



Option



Option



MODEL ECB

- Ideally suited for mobile use
- Battery operation (6 x 1.5V AA)
- Mains Adapter (Optional)
- Rechargeable battery pack optional, at extra cost
- Weighing plate size 318 x 258 mm

MODEL EOB

- Wall mounting for the display device
- Battery operation (6 x 1.5V AA) possible
- Further weighing unit : lb
- Weighing plate size 310 x 300 mm
- Ideally suited for dispatch department



Cat No.	Model KERN	Read out 'd' g	Weighing range kg	Reproducibility g	Linearity g
1008 2011	ECB 20K 2N	10	20	10	±40
1008 5021	ECB 50K 2N	20	50	20	±80
1009 3510	EOB 35K 10	10	35	10	±20
1009 6020	EOB 60K 20	20	60	20	±40



KERN FCB

Bench Scale



Features

- **Second display** on the rear of the balance
- **PRE-TARE function** for manual subtraction of a known container weight, useful for checking fill-levels
- **UNITS:** Further weighing units see the internet

Technical data

- Backlit LCD display, digit height 25 mm
- Dimensions of weighing plate (stainless steel)
W x D 252 x 228 mm
- Overall dimensions W x D x H — 270 x 345 x 106 mm

Cat No.	Model KERN	Read out 'd' g	Weighing range kg	Reproducibility g	Linearity g	Min. Piece count wt. g/pc
1010 0103	FCB 3 K 0.1	0.1	3	0.1	±0.3	0.2
1010 0106	FCB 6 K 0.5	0.5	6	0.5	±1.5	1.0
1010 0112	FCB 12 K 1	1.0	12	1.0	±3.0	2.0