



**KERN KB**  
Precision balance

- Optional Extra :
- Statistics printer
  - Standard printer
  - Rechargeable battery pack




Ring Shaped draft shield  
only with KB 120-3  
& KB 200-3

Cat No.	Model KERN	Read out 'd' g	Weighing range g	Reproducibility g	Linearity g	Weighing Plate mm
1022 0120	KB 120-3N	0.001	121	0.001	± 0.003	Ø 81
1022 0240	KB 240-3N	0.001	200	0.001	± 0.003	Ø 81
1022 1200	KB 1200-2N	0.01	1210	0.01	± 0.03	130 x 130
1022 2000	KB 2000-2N	0.01	2010	0.01	± 0.03	130 x 130
1022 0010	KB 10000-1N	0.1	10100	0.1	± 0.3	150 x 170



**KERN 572**  
Precision balance

- Optional Extra :
- Statistics printer
  - Standard printer
  - DAkks certificate
  - Rechargeable Battery Pack



572-30,  
572-31 & 572-32  
with Glass Draft  
Shield (optional)

Ring Shaped draft shield  
(standard) only with  
572-30, 572-31 & 572-32

Cat No.	Model KERN	Read out 'd' g	Weighing range g	Reproducibility g	Linearity g	Weighing Plate	
1023 0241	572-30	0.001	240	0.001	± 0.003	Ø 105	
1023 0301	572-31	0.001	300	0.002	± 0.003	Ø 105	
1023 0421	572-32	0.001	420	0.002	± 0.005	Ø 105	
1023 1610	572-33	0.01	1600	0.01	± 0.03	Ø 150	
1023 2410	572-35	0.01	2400	0.01	± 0.03	Ø 150	
1023 3010	572-37	0.01	3000	0.02	± 0.05	Ø 150	
1023 1210	572-45	0.10	12000	0.05	± 0.15	160 x 200	
1023 0024	572-57	0.10	24000	0.10	± 0.30	160 x 200	
1023 0501	572-A05	Glass draft shield for Models 572 with d=0.001 g					