



KERN CH • HCB

Hanging scale



KERN CH

Features :

- Data hold function
- Stainless steel hook
- Tape measure extractable (100 cm)
- Block batteries (9 V)
Auto off after 3 minutes without change of load



KERN HCB

Features :

- Data hold function
- Stainless steel hook
- Stainless steel snap link
- Batteries (3 x 1.5 V AA)
Auto off after 3 minutes without change of load

Cat No.	Model KERN	Read out 'd' g	Weighing range kg
1017 0015	CH 15K 20	20	15
1017 0050	CH 50K 50	50	50
1017 0051	CH 50K 100	100	50
1018 0200	HCB 200K 500	500	200



KERN PFB

Precision balance

Low-cost model based on a concept of easy use

Features

- **Concept of easy use:** All primary functions have their own key on the keypad.
- **Compact size,** practical for small spaces.
- **Capacity display:** A bar lights up to show how much of the weighing range is still available.
- **Level indicator and foot screws** to level the balance precisely, fitted as standard, to give the most accurate weighing results,
- **Draft shield** standard (for models with weighing plate sizes A and B).
Removable metal cover with pipette opening.
Weighing space WxDxH 158x143x64 mm

Technical data

- Backlit LCD display, digit height 15 mm
- Dimensions of weighing plate (stainless steel)

A — Ø 80 mm

B — Ø 120 mm



Cat No.	Model KERN	Read out 'd' g	Weighing range g	Reproducibility g	Linearity g	Weighing Plate mm
1029 0120	PFB 120-3	0.001	120	0.001	± 0.003	Ø 80
1029 0200	PFB 200-3	0.001	200	0.001	± 0.003	Ø 80
1029 0300	PFB 300-3	0.001	300	0.001	± 0.003	Ø 80
1029 1200	PFB 1200-2	0.01	1200	0.01	± 0.03	Ø 120

