



Test weights :

Single as well as set of test weights are available in stainless steel or nickel plated brass. They can be supplied with German DAkkS certificate or locally with ERTL traceability, at extra cost.

Single weight Stainless Steel Class, E2

Cat. No.	Description
1000 0535	1 g in box
1000 0536	2 g in box
1000 0537	5 g in box
1000 0538	10 g in box
1000 0539	20 g in box
1000 0540	50 g in box
1000 0541	100 g in box
1000 0542	200 g in box
1000 0543	500 g in box
1000 0544	1 kg in box
1000 0545	2 kg in box
1000 0546	5 kg in box
1000 0547	10 kg in box
1000 0548	20 kg in box
1000 0549	50 kg in box

Single weight Stainless Steel, Class F1

Cat. No.	Description
1000 0551	1 g in box
1000 0552	2 g in box
1000 0553	5 g in box
1000 0554	10 g in box
1000 0555	20 g in box
1000 0556	50 g in box
1000 0557	100 g in box
1000 0558	200 g in box
1000 0559	500 g in box
1000 0560	1 kg in box
1000 0561	2 kg in box
1000 0562	5 kg in box
1000 0563	10 kg in box
1000 0564	20 kg in box
1000 0565	50 kg in box

Weight Set, Stainless Steel, in plastic carrying case, Class E2

Cat. No.	Description
1000 0571	Range 1 mg - 200 g

Weight Set, Stainless Steel, in plastic carrying case, Class F1

Cat. No.	Description
1000 0581	Range 1 mg - 200 g

Note : Weight boxes and weights in other ranges and class E1, E2, F1, F2, M1, M2 & M3 are also available. Details / Price on Request