


pH Tutor		
Measuring Parameter	pH/°C	
Highlights	Dual display standard bench	
pH	Range	0.00 to 14.00 pH
	Resolution	0.01pH
	Accuracy	±0.01pH
	Calibration Points	Up to 3
	Buffer Sets	USA, NIST
Temperature	Range	0 to 100°C
	Resolution	0.1 °C
	Accuracy	±0.3 °C
Meter Features	Temperature Compensation	ATC/MTC (0 to 100°C)
	Input	BNC, phono
	Power Requirements	9 V DC Adapter, 500 mA (220 VAC)
Dimensions (LxWxH); Weight	Meter	23 x 18 x 6 cm ; 570 g
	Boxed	40 x 26 x 9 cm ; 1200 g

Order Code	ECPHTUTOR-DS	ECPHTUTOR-D	ECPHTUTOR-S	ECPHTUTOR-I
pH Electrode	EUTECH "3-in-1" pH/Temp Order Code ECFC7352901B	EUTECH "3-in-1" pH/Temp Order Code ECFC7352901B	EUTECH Electrode Order Code ECFC7252101B	Indigenous pH Electrode Order Code AN-786-1
Stand	EUTECH electrode stand with swivel arm	Indigenous electrode stand	Indigenous electrode stand	Indigenous electrode stand



Product Code	Cat No.	Description
ECPHTUTOR-DS	8034 011	CyberScan pH Tutor meter with plastic-body "3-in-1" pH/Temp.(ATC in-built) combination electrode, 12x90mm, BNC connector, 1m cable length, integrated electrode stand with swivel arm and power adapter
ECPHTUTOR-D	8034 001	CyberScan pH tutor pH meter with plastic-body "3-in-1" pH/Temp.(ATC in-built) ECFC7352901 combination electrode, 12x90mm, BNC connector, 1m cable length, indigenous electrode holder and power adapter
ECPHTUTOR-S	8034 002	CyberScan pH Tutor pH meter with single junction pH electrode ECFC7252101B combination electrode, 12x90mm, BNC connector, 1 m cable length, indigenous electrode holder and power adapter. ATC probe to be ordered separately
ECPHTUTOR-I	8034 010	CyberScan pH Tutor meter with indigenous refillable glass pH electrode AN-786-1, BNC connector, 1m cable length, indigenous electrode stand and power adapter. ATC probe to be ordered separately
ECFC7252101B	8034 004	General purpose plastic-body single junction gel-Filled pH combination electrode, 12x90mm, BNC connector, 1m cable length
ECFG7350401B	8034 006	General purpose glass-body single junction refillable pH combination electrode, annular ceramic reference junction with protective sensor guard, 12x110mm, BNC connector, 1m cable length & 10ml refilling electrolyte
ECFC7352901B	8034 005	General purpose plastic-body "3-in-1" pH/Temp. combination electrode, 12x90mm, BNC connector, 1m cable length (for pH Tutor)
AN-786-1	8034 009	General purpose glass-body refillable indigenous pH combination electrode, BNC connector, 1m cable length, Make ANTECH.
ECPH5TEMB01P	8034 003	ATC Probe, 1m cable length (for pH Tutor)
ECELSTDCNEW	8034 007	Electrode stand with swivel arm (for pH Tutor)
60X030130	8050 007	100 /240VAC SMPS Power Adapter with US/UK/EUR/Japan plug, 9V,6W, CENTRE + ve (for Tutor)